


		Genotyp	ing Report		· · < *
Strain ID	T005478	Strain Type	CKO(Cas9)	Genetic Background	C57BL/6JGpt
Designer	Ya'nan Xu	Gene Name	~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~	Mapk12	Co l
Strategy of (Genotyping	24		armar.	~< x




Note: The sizes of WT and Targeted band are shown below.

2. Primer Information

PCR No.	Primer No.	Sequence	Band Size	
(1)(5'arm)	T005478-F1	CTCAATGGAGCAAAGCCAAAG	WT:316bp	
	T005478-R1	TGAGGGAGTTGGACAACAGATAGAC	Targeted:373bp	
2)(3'arm)	T005478-F2	TAGCATCCTCATGTCCAGACCAC	WT:291bp	
	T005478-R2 TAGGCAAACTGAGAGAGTTCTCTGAG		- Targeted:348bp	

3. Gel Image & Conclusion





江 就 集 萃 药 康 生 物 科 技 股 份 有 限 公 司 Gem Pharmatech Co.,Ltd

product band position and size meet the theoretical requirements.

⁽²⁾ Control (B) : PCR amplification was performed without template in the PCR reagent to monitor whether the reagent was contaminated.

4. PCR Condition

PCR Reaction	Component		M. M		
Seg.	reaction co	reaction component			
1^{γ}	2 × Rapid Taq Master Mix (Vazym	2 × Rapid Taq Master Mix (Vazyme P222)			
2 73	ddH2O		9.5		
3	Primer A(10pmol/µl)	. 3/x	1		
1	Primer B(10pmol/μl)	Primer B(10pmol/µl)			
5	Template(≈100ng/µl)				
PCR program	① priority selection	í Constantino de Cons			
Seg.	Temp.	Time	Cycle		
1	95℃	5min	CULE .		
2 6	98°C	30s	20×		
3 ^N S,	65℃*(-0.5℃/cycle)	30s			
4 7	72°C	45s*			
5	98°C	30s	20×		
5 6	55℃*	30s			
	72°C	45s*			
3	72°C	5min	Dr. Slx		
Ð	10 ℃	hold	() ₂ , (9		
PCR program	② the second choice	Tax Con	2 ² C/		
Seg.	Temp.	Time	Cycle		
1 97	95℃	5min	AND - CA		
2	98°C	30s	35×		
3	58℃*	58℃* 30s			
4	72°C	45s*			
5 7/	72℃	5min			
6	10°C	hold	ngr.		

Note*: Annealing temperature and extension time can be determined according to the actual amplification situation and amplification enzyme efficiency.