


Genotyping Report

Strain ID	T027634	Strain Type	KO(Cas9)	Genetic Background	C57BL/6JGpt
Designer	Zifan Lin	Gene Name	<i>Kmt5a</i>		

1. Strategy of Genotyping




2. Primer Information

PCR No.	Primer No.	Primer Name	Sequence	Band Size
PCR①	F1	JS08567-Kmt5a-5wt-tF1	GTTTGACTGTGGCCTTGTGGC	WT:2262 bp KO:266bp
	R1	JS08567-Kmt5a-3wt-tR1	CAACTGACTGATGTGTCACCCAGG	
PCR②	F2	JS18567-Kmt5a-wt-F1	GCCAGACTTATCCCTTAGGGACTC	WT: 314bp KO:0 bp
	R2	JS18567-Kmt5a-wt-R1	GCGCATTGCAGAGCACTTGTT	

3. Gel Image

ccggagttacttggtgtctagtggcttacc---1996bp---gaattgaggagctatgttgtctgagccaca
--



Note: P: Positive control; W: Wildtype control; B: Blank control (ddH₂O); M: DNA Ladder

4. PCR Condition

PCR Reaction Component			
Seg.	reaction component		Volume (μl)
1	2 × Taq Master Mix ， Dye Plus， (Vazyme P112-03)		12.5
2	ddH2O		9.5
3	Primer A(10pmol/μl)		1
4	Primer B(10pmol/μl)		1
5	Template(≈100ng/μl)		1
PCR program			
Seg.	Temp.	Time	Cycle
1	95℃	5min	20×
2	98℃	30s	
3	65℃ （-0.5℃/cycle）	30s	

4	72℃	45s	20×
5	98℃	30s	
6	55℃	30s	
7	72℃	45s	
8	72℃	5min	
9	10℃	hold	

Reviewer:Ting Sun

Date:2020-7-1