


Genotyping Report

Strain ID	T013798	Strain Type	KO(Cas9)	Genetic Background	C57BL/6JGpt
Designer	Zifan Lin	Gene Name	Mul1		

1. Strategy of Genotyping




2. Primer Information

PCR No.	Primer No.	Primer Name	Sequence	Band Size
PCR①	F1	JS04252-Mul1-5wt-tF1	AACCCACACCCCTTAGTGACACAC	WT:6001bp KO:275bp
	R1	JS04252-Mul1-3wt-tR1	ACAGGATGAAGAGGTGGCAAC	
PCR②	F2	JS14252-Mul1-wt-tF1	TTTGCTGGGAGCACTTGAGAC	WT:365bp KO:0bp
	R2	JS14252-Mul1-wt-tR1	AGAGAGGGCTCTTTCACTCCATG	

3. Gel Image

tcatttatccactcagaaaattccatgag-----5726bp-----agagggtgtgggcatttgaataaccagctta
--



4. PCR Condition

PCR Reaction Component			
Seg.	reaction component		Volume (μl)
1	2 × Taq Master Mix , Dye Plus, (Vazyme P112-03)		12.5
2	ddH2O		9.5
3	Primer A(10pmol/μl)		1
4	Primer B(10pmol/μl)		1
5	Template(≈100ng/μl)		1
PCR program			
Seg.	Temp.	Time	Cycle
1	95℃	5min	20×
2	98℃	30s	
3	65℃ （-0.5℃/cycle）	30s	
4	72℃	45s	
5	98℃	30s	20×
6	55℃	30s	
7	72℃	45s	
8	72℃	5min	
9	10℃	hold	

Reviewer:Ting Sun

Date:2020-7-1