


## Genotyping Report

|           |          |             |                |                    |             |
|-----------|----------|-------------|----------------|--------------------|-------------|
| Strain ID | T028473  | Strain Type | KO(Cas9)       | Genetic Background | C57BL/6JGpt |
| Designer  | Yin Chen | Gene Name   | <i>B4galt1</i> |                    |             |

### 1. Strategy of Genotyping




Wild type: ①PCR reaction obtains a single WT band; ②PCR reaction obtains a single WT band.  
Heterozygote: ①PCR reaction obtains a WT band and a KO band; ②PCR reaction obtains a WT band.  
Homozygote: ①PCR reaction obtains a single KO band; ②PCR reaction without product.  
Note: 1)The sizes of WT and Targeted band are shown below.  
2)If the WT band is too large, it may not be possible to obtain a WT band.

### 2. Primer Information

| PCR No. | Primer No. | Primer Name             | Sequence                        | Band Size              |
|---------|------------|-------------------------|---------------------------------|------------------------|
| PCR①    | F1         | JS09406-B4galt1-5wt-tF1 | AAAGCTGGCCACCGCAGAATAAC         | WT:3250bp<br>KO: 370bp |
|         | R1         | JS09406-B4galt1-3wt-tR1 | CAACCCAGCGTTATTCATGTGACC        |                        |
| PCR②    | F2         | JS19406-B4galt1-wt-F1   | ACTATGGCATCTACGTCATCAATCA<br>GG | WT:240bp<br>KO:0bp     |
|         | R2         | JS19406-B4galt1-wt-R1   | AGCACACCGTATCTATTTCTCTTGG<br>C  |                        |

### 3. Gel Image

|                                                                                              |
|----------------------------------------------------------------------------------------------|
| gaggggcttcttggccgcacagtttttcagtgcc----- <b>2880bp</b> -----cttagtatcagaaaaggactcttttaggagctg |
|----------------------------------------------------------------------------------------------|



Note: P: Positive control; WT: Wildtype control; B: Blank control (ddH<sub>2</sub>O); M: DNA Ladder

#### 4. PCR Condition

| PCR Reaction Component |                                                 |      |             |
|------------------------|-------------------------------------------------|------|-------------|
| Seg.                   | reaction component                              |      | Volume (μl) |
| 1                      | 2 × Taq Master Mix , Dye Plus, (Vazyme P112-03) |      | 12.5        |
| 2                      | ddH2O                                           |      | 9.5         |
| 3                      | Primer A(10pmol/μl)                             |      | 1           |
| 4                      | Primer B(10pmol/μl)                             |      | 1           |
| 5                      | Template(≈100ng/μl)                             |      | 1           |
| PCR program            |                                                 |      |             |
| Seg.                   | Temp.                                           | Time | Cycle       |
| 1                      | 95℃                                             | 5min | 20×         |
| 2                      | 98℃                                             | 30s  |             |
| 3                      | 65℃ （-0.5℃/cycle）                               | 30s  |             |
| 4                      | 72℃                                             | 45s  |             |
| 5                      | 98℃                                             | 30s  | 20×         |
| 6                      | 55℃                                             | 30s  |             |
| 7                      | 72℃                                             | 45s  |             |
| 8                      | 72℃                                             | 5min |             |
| 9                      | 10℃                                             | hold |             |

Reviewer:Ting Sun

Date:2020-7-1