


Genotyping Report

Strain ID	T017886	Strain Type	KO(Cas9)	Genetic Background	C57BL/6JGpt
Designer	Zifan Lin	Gene Name	Ift80		

1. Strategy of Genotyping




2. Primer Information

PCR No.	Primer No.	Primer Name	Sequence	Band Size
PCR①	F1	JS06099-Ift80-5wt-tF1	CACTTGAAAGTGTTGGGTTTTTCAGT TG	WT:25136bp KO: 236bp
	R1	JS06099-Ift80-3wt-tR1	GCTCCACTGTCATTATCAAAACATG AGA	
PCR②	F2	JS16099-Ift80-wt-tF1	AGGATCTTGCAGAGTGACCCGCTT	WT:265bp KO:0bp
	R2	JS16099-Ift80-wt-tR1	ATGTGATACAGTGGCTGTCTTTGCC	

3. Gel Image

TTGAAAAGTCTGCCTTTTATGTGGC-----24900bp-----CAAGCAAGTCATGAGAGACAAGCCA



4. PCR Condition

PCR Reaction Component			
Seg.	reaction component		Volume (μl)
1	2 × Taq Master Mix , Dye Plus, (Vazyme P112-03)		12.5
2	ddH2O		9.5
3	Primer A(10pmol/μl)		1
4	Primer B(10pmol/μl)		1
5	Template(≈100ng/μl)		1
PCR program			
Seg.	Temp.	Time	Cycle
1	95℃	5min	
2	98℃	30s	20×

3	65℃ (-0.5℃/cycle)	30s	20×
4	72℃	45s	
5	98℃	30s	
6	55℃	30s	
7	72℃	45s	
8	72℃	5min	
9	10℃	hold	

Reviewer:Ting Sun

Date:2020-7-1