


## Genotyping Report

|           |           |             |          |                    |             |
|-----------|-----------|-------------|----------|--------------------|-------------|
| Strain ID | T014504   | Strain Type | KO(Cas9) | Genetic Background | C57BL/6JGpt |
| Designer  | Zifan Lin | Gene Name   | Mettl14  |                    |             |

### 1. Strategy of Genotyping



### 2. Primer Information

| PCR No. | Primer No. | Primer Name             | Sequence                  | Band Size               |
|---------|------------|-------------------------|---------------------------|-------------------------|
| PCR①    | F1         | JS03443-Mettl14-5wt-tF1 | GCCGTGTTTGTACCTAATACTGC   | WT:13502bp<br>KO: 317bp |
|         | R1         | JS03443-Mettl14-3wt-tR1 | ACAGCAATGCACAGCTTAGCTG    |                         |
| PCR②    | F2         | JS13443-Mettl14-wt-tF1  | GCCTATTGATAGTGTGCTAAGGGTG | WT:379bp<br>KO:0bp      |
|         | R2         | JS13443-Mettl14-wt-tR1  | CAGTTCTCAAGAGGTTTCAGGCAG  |                         |

### 3. Gel Image

|                                                                                      |
|--------------------------------------------------------------------------------------|
| attcgacagaattcctccagcgccactta----- <b>13185bp</b> -----ttttggcacctcctaattatccgaatttt |
|--------------------------------------------------------------------------------------|



Note: P: Positive control; WT: Wildtype control; B: Blank control (ddH<sub>2</sub>O); M: DNA Ladder

#### 4. PCR Condition

| PCR Reaction Component |                                                 |      |             |
|------------------------|-------------------------------------------------|------|-------------|
| Seg.                   | reaction component                              |      | Volume (μl) |
| 1                      | 2 × Taq Master Mix , Dye Plus, (Vazyme P112-03) |      | 12.5        |
| 2                      | ddH2O                                           |      | 9.5         |
| 3                      | Primer A(10pmol/μl)                             |      | 1           |
| 4                      | Primer B(10pmol/μl)                             |      | 1           |
| 5                      | Template(≈100ng/μl)                             |      | 1           |
| PCR program            |                                                 |      |             |
| Seg.                   | Temp.                                           | Time | Cycle       |
| 1                      | 95℃                                             | 5min | 20×         |
| 2                      | 98℃                                             | 30s  |             |
| 3                      | 65℃ （-0.5℃/cycle）                               | 30s  |             |
| 4                      | 72℃                                             | 45s  |             |
| 5                      | 98℃                                             | 30s  | 20×         |
| 6                      | 55℃                                             | 30s  |             |
| 7                      | 72℃                                             | 45s  |             |
| 8                      | 72℃                                             | 5min |             |
| 9                      | 10℃                                             | hold |             |

Reviewer:Ting Sun

Date:2020-7-1