


## Genotyping Report

|           |           |             |              |                    |             |
|-----------|-----------|-------------|--------------|--------------------|-------------|
| Strain ID | T012282   | Strain Type | KO(Cas9)     | Genetic Background | C57BL/6JGpt |
| Designer  | Zifan Lin | Gene Name   | <i>Pold3</i> |                    |             |

### 1. Strategy of Genotyping





### 2. Primer Information

| PCR No. | Primer No. | Sequence                  | Band Size               |
|---------|------------|---------------------------|-------------------------|
| PCR①    | T012282-F1 | GTGTCAACCACCACACTGTGAAG   | WT:12296bp<br>KO: 304bp |
|         | T012282-R1 | GATTACAAACACATGACACAACGAC |                         |
| PCR②    | T012282-F2 | GAGAAGGCTGAGTAGGTTTCCATC  | WT:360bp<br>KO:0bp      |
|         | T012282-R2 | CTGACCAAGTTCAAGGAGCAGTAAC |                         |

### 3. Gel Image

|                                                                                          |
|------------------------------------------------------------------------------------------|
| cagggtttgtgatgctgtgctgactgggcccctca-----11994bp-----ACgtgcataggagggtgctgctgctggaacacgtca |
|------------------------------------------------------------------------------------------|



Note: P: Positive control; WT: Wildtype control; B: Blank control (ddH<sub>2</sub>O); M: DNA Ladder

#### 4. PCR Condition

| PCR Reaction Component |                                                 |      |             |
|------------------------|-------------------------------------------------|------|-------------|
| Seg.                   | reaction component                              |      | Volume (μl) |
| 1                      | 2 × Taq Master Mix , Dye Plus, (Vazyme P112-03) |      | 12.5        |
| 2                      | ddH2O                                           |      | 9.5         |
| 3                      | Primer A(10pmol/μl)                             |      | 1           |
| 4                      | Primer B(10pmol/μl)                             |      | 1           |
| 5                      | Template(≈100ng/μl)                             |      | 1           |
| PCR program            |                                                 |      |             |
| Seg.                   | Temp.                                           | Time | Cycle       |
| 1                      | 95℃                                             | 5min | 20×         |
| 2                      | 98℃                                             | 30s  |             |
| 3                      | 65℃ (-0.5℃/cycle)                               | 30s  |             |
| 4                      | 72℃                                             | 45s  |             |
| 5                      | 98℃                                             | 30s  | 20×         |
| 6                      | 55℃                                             | 30s  |             |
| 7                      | 72℃                                             | 45s  |             |
| 8                      | 72℃                                             | 5min |             |
| 9                      | 10℃                                             | hold |             |

Reviewer:Ting Sun

Date:2021-7-22