


Genotyping Report

Strain ID	T014843	Strain Type	KO(Cas9)	Genetic Background	C57BL/6JGpt
Designer	Zifan Lin	Gene Name	Foxk1		

1. Strategy of Genotyping




2. Primer Information

PCR No.	Primer No.	Sequence	Band Size
PCR①	T014843-F1	ACTCGTCCTGTGCCTTACCC	WT:15334bp KO: 312bp
	T014843-R1	CTAACTATTGAGGAGGTGGAGGC	
PCR②	T014843-F2	ATGTAGTGAGTTCGAGGCCAAGC	WT:352bp KO:0bp
	T014843-R2	ACAGCACAACCTCTGGCTTTGCTG	

3. Gel Image

aggcagcctagctcagcatgagcctctggccccc----- 15022bp -----tgatggfactgggaatcagaccaggagtgagtct
--



4. PCR Condition

PCR Reaction Component			
Seg.	reaction component		Volume (μl)
1	2 × Taq Master Mix , Dye Plus, (Vazyme P112-03)		12.5
2	ddH2O		9.5
3	Primer A(10pmol/μl)		1
4	Primer B(10pmol/μl)		1
5	Template(≈100ng/μl)		1
PCR program			
Seg.	Temp.	Time	Cycle
1	95℃	5min	20×
2	98℃	30s	
3	65℃ （-0.5℃/cycle）	30s	
4	72℃	45s	
5	98℃	30s	20×
6	55℃	30s	
7	72℃	45s	
8	72℃	5min	
9	10℃	hold	

Reviewer:Ting Sun

Date:2021-7-26