


Genotyping Report

Strain ID	T014846	Strain Type	KO(Cas9)	Genetic Background	C57BL/6JGpt
Designer	Zifan Lin	Gene Name	Psmc2		

1. Strategy of Genotyping



2. Primer Information

PCR No.	Primer No.	Sequence	Band Size
PCR①	T014846-F1	ACCCTCAGTTGAGTCACATTTGTC	WT:2884bp KO: 226bp
	T014846-R1	AGCCCTAGCTGTACTAGAACTTGC	
PCR②	T014846-F2	AGAGACGTAGCAGTTTCTGATGCG	WT:262bp KO:0bp
	T014846-R2	AACCCTTCAGGACCAAGAGGAG	

3. Gel Image

tcctaattggccttgctttgtccttgactcccaca----- 2658bp -----ggtagagggcacaccctttgatccaacacttgg



4. PCR Condition

PCR Reaction Component			
Seg.	reaction component		Volume (μl)
1	2 × Taq Master Mix , Dye Plus, (Vazyme P112-03)		12.5
2	ddH2O		9.5
3	Primer A(10pmol/μl)		1
4	Primer B(10pmol/μl)		1
5	Template(≈100ng/μl)		1
PCR program			
Seg.	Temp.	Time	Cycle
1	95℃	5min	20×
2	98℃	30s	
3	65℃ (-0.5℃/cycle)	30s	
4	72℃	45s	
5	98℃	30s	20×
6	55℃	30s	
7	72℃	45s	
8	72℃	5min	
9	10℃	hold	



Reviewer:Ting Sun

Date:2021-7-26