


Genotyping Report

Strain ID	T014290	Strain Type	KO(Cas9)	Genetic Background	C57BL/6JGpt
Designer	Zifan Lin	Gene Name	<i>Rnf13</i>		

1. Strategy of Genotyping




2. Primer Information

PCR No.	Primer No.	Sequence	Band Size
PCR①	T014290-F1	GCTCCGTGTTTGATACTTCTGAC	WT:29663bp KO: 402bp
	T014290-R1	TTCCCTGTGCTGCATCGTAC	
PCR②	T014290-F2	AATCTTAGCACTCAGGAGGATGG	WT:353bp KO:0bp
	T014290-R2	GCTGGGATATGCAATAGTGTCTG	

3. Gel Image

aaattactttatgacttgaggagtcctacag--- 29261bp ---ctcgcatatgctttcttatcataggtggcc



Note: P: Positive control; WT: Wildtype control; B: Blank control (ddH₂O); M: DNA Ladder

4. PCR Condition

PCR Reaction Component			
Seg.	reaction component		Volume (μl)
1	2 × Taq Master Mix , Dye Plus, (Vazyme P112-03)		12.5
2	ddH2O		9.5
3	Primer A(10pmol/μl)		1
4	Primer B(10pmol/μl)		1
5	Template(≈100ng/μl)		1
PCR program			
Seg.	Temp.	Time	Cycle
1	95℃	5min	20×
2	98℃	30s	
3	65℃ (-0.5℃/cycle)	30s	
4	72℃	45s	
5	98℃	30s	20×
6	55℃	30s	
7	72℃	45s	
8	72℃	5min	
9	10℃	hold	

Reviewer:Ting Sun

Date:2021-7-26