


Genotyping Report

Strain ID	T013778	Strain Type	KO(Cas9)	Genetic Background	C57BL/6JGpt
Designer	Zifan Lin	Gene Name	<i>Kcnn2</i>		

1. Strategy of Genotyping



2. Primer Information

PCR No.	Primer No.	Sequence	Band Size
PCR①	T013778-F1	TTAACCCTGTCTGCGCCGTA	WT:3614bp KO: 236bp
	T013778-R1	AGGAAGAAGTCTCAGTGGTGCAGG	
PCR②	T013778-F2	AAGCATGTTAAGTGGCTCCTGGAG	WT:343bp KO:0bp
	T013778-R2	AGCACTTCCAATCTCTTTCCAGGC	

3. Gel Image

ggacagatccctgcctgggtcggatgtgatctct----- 3378bp -----caggtcgggaaggggccagcatgaaagtaggcact
--



4. PCR Condition

PCR Reaction Component			
Seg.	reaction component		Volume (μl)
1	2 × Taq Master Mix , Dye Plus, (Vazyme P112-03)		12.5
2	ddH2O		9.5
3	Primer A(10pmol/μl)		1
4	Primer B(10pmol/μl)		1
5	Template(≈100ng/μl)		1
PCR program			
Seg.	Temp.	Time	Cycle
1	95℃	5min	20×
2	98℃	30s	
3	65℃ （-0.5℃/cycle）	30s	
4	72℃	45s	
5	98℃	30s	20×
6	55℃	30s	
7	72℃	45s	
8	72℃	5min	
9	10℃	hold	

Reviewer:Ting Sun

Date:2021-7-26



江苏集萃药康生物科技股份有限公司
GemPharmatech Co.,Ltd