


Genotyping Report

Strain ID	T027968	Strain Type	KO(Cas9)	Genetic Background	C57BL/6JGpt
Designer	Zifan Lin	Gene Name	<i>Pelp1</i>		

1. Strategy of Genotyping




2. Primer Information

PCR No.	Primer No.	Sequence	Band Size
PCR①	T027968-F1	GCTATTACCAGCACCTGACTTCAC	WT:1775bp KO:276bp
	T027968-R1	CTCAGAGAGAGGTAGAGTAGCAAGC	
PCR②	T027968-F2	CCTTTGGGAGTCGATTATCCTAG	WT:451bp KO:0bp
	T027968-R2	TCACCTCAGCAGACCTTCCTTAG	

3. Gel Image

tcgtgtcttggttaggagcatttgggttgctag----- 1499bp -----ccatggattcaagagatttttccattccaggag



4. PCR Condition

PCR Reaction Component			
Seg.	reaction component		Volume (μl)
1	2 × Taq Master Mix , Dye Plus, (Vazyme P112-03)		12.5
2	ddH2O		9.5
3	Primer A(10pmol/μl)		1
4	Primer B(10pmol/μl)		1
5	Template(≈100ng/μl)		1
PCR program			
Seg.	Temp.	Time	Cycle
1	95℃	5min	20×
2	98℃	30s	
3	65℃ （-0.5℃/cycle）	30s	
4	72℃	45s	
5	98℃	30s	20×
6	55℃	30s	
7	72℃	45s	
8	72℃	5min	
9	10℃	hold	

Reviewer:Ting Sun

Date:2021-7-10