

***Zfp853* Cas9-CKO Strategy**

Designer: Yumeng Wang

Reviewer: Longyun Hu

Design Date: 2020-10-21

Project Overview

Project Name

Zfp853

Project type


Cas9-CKO

Strain background

C57BL/6JGpt

Conditional Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Zfp853* gene. The schematic diagram is as follows:



Technical routes

- The *Zfp853* gene has 3 transcripts. According to the structure of *Zfp853* gene, exon1-exon3 of *Zfp853*-202(ENSMUST00000212355.1) transcript is recommended as the knockout region. The region contains all of the coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Zfp853* gene. The brief process is as follows: CRISPR/Cas9 system and Donor were microinjected into the fertilized eggs of C57BL/6JGpt mice. Fertilized eggs were transplanted to obtain positive F0 mice which were confirmed by PCR and sequencing. A stable F1 generation mouse model was obtained by mating positive F0 generation mice with C57BL/6JGpt mice.
- The flox mice will be knocked out after mating with mice expressing Cre recombinase, resulting in the loss of function of the target gene in specific tissues and cell types.

Notice

- The *Zfp853* gene is located on the Chr5. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of loxp insertion on gene transcription, RNA splicing and protein translation cannot be predicted at existing technological level.

Gene information (NCBI)

Zfp853 zinc finger protein 853 [Mus musculus (house mouse)]

Gene ID: 330230, updated on 13-Mar-2020

Summary




Official Symbol	Zfp853 provided by MGI
Official Full Name	zinc finger protein 853 provided by MGI
Primary source	MGI:MGI:2685638
See related	Ensembl:ENSMUSG00000093910
Gene type	protein coding
RefSeq status	VALIDATED
Organism	Mus musculus
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	Gm792
Expression	Broad expression in ovary adult (RPKM 1.8), limb E14.5 (RPKM 1.7) and 16 other tissues See more
Orthologs	human all

Transcript information (Ensembl)


The gene has 3 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Zfp853-203	ENSMUST00000212715.1	3133	718aa	Protein coding	-	A0A1D5RMN1	TSL:5 GENCODE basic APPRIS ALT2
Zfp853-202	ENSMUST00000212355.1	2871	733aa	Protein coding	-	A0A1D5RM95	TSL:5 GENCODE basic APPRIS P5
Zfp853-201	ENSMUST00000180336.1	1956	652aa	Protein coding	-	J3QPU9	TSL:5 GENCODE basic APPRIS ALT2

The strategy is based on the design of *Zfp853-202* transcript,the transcription is shown below:



Genomic location distribution



集萃药康
GemPharmatech



If you have any questions, you are welcome to inquire.

Tel: 400-9660890

