

Kif24 Cas9-KO Strategy

Designer: Miaomiao Cui

Reviewer: Lingyan Wu

Design Date: 2021-5-31

Project Overview

Project Name

Kif24

Project type


Cas9-KO

Strain background

C57BL/6JGpt

Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Kif24* gene. The schematic diagram is as follows:



- The *Kif24* gene has 8 transcripts. According to the structure of *Kif24* gene, exon2 of *Kif24-202*(ENSMUST00000108055.9) transcript is recommended as the knockout region. The region contains start codon ATG. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Kif24* gene. The brief process is as follows: sgRNA was transcribed in vitro. Cas9 and sgRNA were microinjected into the fertilized eggs of C57BL/6JGpt mice. Fertilized eggs were transplanted to obtain positive F0 mice which were confirmed by PCR and sequencing. A stable F1 generation mouse model was obtained by mating positive F0 generation mice with C57BL/6JGpt mice.

- The *Kif24* gene is located on the Chr4. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of the gene knockout on gene transcription, RNA splicing and protein translation cannot be predicted at the existing technology level.

Gene information (NCBI)

Kif24 kinesin family member 24 [Mus musculus (house mouse)]

Gene ID: 109242, updated on 17-Feb-2021

Summary




Official Symbol	Kif24 provided by MGI
Official Full Name	kinesin family member 24 provided by MGI
Primary source	MGI:MGI:1918345
See related	Ensembl:ENSMUSG00000028438
Gene type	protein coding
RefSeq status	VALIDATED
Organism	Mus musculus
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	4933425J19Rik, 9430029L23Rik
Expression	Broad expression in testis adult (RPKM 3.9), CNS E11.5 (RPKM 2.1) and 19 other tissues See more
Orthologs	human all

Transcript information (Ensembl)


The gene has 8 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Kif24-202	ENSMUST00000108055.9	5593	1356aa	Protein coding	CCDS18060		TSL:1 , GENCODE basic , APPRIS P1 ,
Kif24-201	ENSMUST00000030148.6	4333	1222aa	Protein coding	-		TSL:1 , GENCODE basic ,
Kif24-208	ENSMUST00000154535.8	731	141aa	Protein coding	-		CDS 3' incomplete , TSL:3 ,
Kif24-207	ENSMUST00000154175.8	1756	No protein	Processed transcript	-		TSL:1 ,
Kif24-204	ENSMUST00000131532.2	1151	No protein	Processed transcript	-		TSL:1 ,
Kif24-206	ENSMUST00000141951.2	754	No protein	Processed transcript	-		TSL:3 ,
Kif24-203	ENSMUST00000127376.2	705	No protein	Processed transcript	-		TSL:3 ,
Kif24-205	ENSMUST00000133497.2	623	No protein	Processed transcript	-		TSL:5 ,


The strategy is based on the design of *Kif24-202* transcript,the transcription is shown below:



Genomic location distribution



Protein domain



If you have any questions, you are welcome to inquire.

Tel: 025-5864 1534

