

Slc35b4 Cas9-CKO Strategy

Designer: Huan Wang

Reviewer: Shanhong Tao

Design Date: 2021-3-23

Project Overview

Project Name

Slc35b4

Project type


Cas9-CKO

Strain background

C57BL/6JGpt

Conditional Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Slc35b4* gene. The schematic diagram is as follows:



Technical routes

- The *Slc35b4* gene has 4 transcripts. According to the structure of *Slc35b4* gene, exon2-exon3 of *Slc35b4-201*(ENSMUST00000019143.8) transcript is recommended as the knockout region. The region contains 217bp coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Slc35b4* gene. The brief process is as follows: CRISPR/Cas9 system and Donor were microinjected into the fertilized eggs of C57BL/6JGpt mice. Fertilized eggs were transplanted to obtain positive F0 mice which were confirmed by PCR and sequencing. A stable F1 generation mouse model was obtained by mating positive F0 generation mice with C57BL/6JGpt mice.
- The flox mice will be knocked out after mating with mice expressing Cre recombinase, resulting in the loss of function of the target gene in specific tissues and cell types.

- The *Slc35b4* gene is located on the Chr6. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of loxp insertion on gene transcription, RNA splicing and protein translation cannot be predicted at existing technological level.

Gene information (NCBI)

Slc35b4 solute carrier family 35, member B4 [Mus musculus (house mouse)]

Gene ID: 58246, updated on 13-Mar-2020

Summary




Official Symbol	Slc35b4 provided by MGI
Official Full Name	solute carrier family 35, member B4 provided by MGI
Primary source	MGI:MGI:1931249
See related	Ensembl:ENSMUSG00000018999
Gene type	protein coding
RefSeq status	PROVISIONAL
Organism	Mus musculus
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	4930474D06Rik, AB041549, AW494345
Expression	Ubiquitous expression in kidney adult (RPKM 32.1), genital fat pad adult (RPKM 17.7) and 27 other tissues See more
Orthologs	human all

Transcript information (Ensembl)


The gene has 4 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Slc35b4-201	ENSMUST00000019143.8	6253	331aa	Protein coding	CCDS19989	Q8CIA5	TSL:1 GENCODE basic APPRIS P1
Slc35b4-202	ENSMUST00000138250.1	984	41aa	Nonsense mediated decay	-	S4R2P0	TSL:5
Slc35b4-203	ENSMUST00000145857.7	4110	No protein	Retained intron	-	-	TSL:2
Slc35b4-204	ENSMUST00000146514.7	2320	No protein	Retained intron	-	-	TSL:1


The strategy is based on the design of *Slc35b4-201* transcript,the transcription is shown below:




Genomic location distribution



Protein domain



Mouse phenotype description(MGI)



Phenotypes affected by the gene are marked in blue. Data quoted from MGI database(<http://www.informatics.jax.org/>).

If you have any questions, you are welcome to inquire.

Tel: 400-9660890

