

Pigs Cas9-CKO Strategy

Designer: Longyun Hu

Reviewer: Rui Xiong

Design Date: 2021-3-29

Project Overview

Project Name

Pigs

Project type


Cas9-CKO

Strain background

C57BL/6JGpt

Conditional Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Pigs* gene. The schematic diagram is as follows:



- The *Pigs* gene has 1 transcript. According to the structure of *Pigs* gene, exon4-exon5 of *Pigs*-201(ENSMUST00000048073.9) transcript is recommended as the knockout region. The region contains 182bp coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Pigs* gene. The brief process is as follows: CRISPR/Cas9 system and Donor were microinjected into the fertilized eggs of C57BL/6JGpt mice. Fertilized eggs were transplanted to obtain positive F0 mice which were confirmed by PCR and sequencing. A stable F1 generation mouse model was obtained by mating positive F0 generation mice with C57BL/6JGpt mice.
- The flox mice will be knocked out after mating with mice expressing Cre recombinase, resulting in the loss of function of the target gene in specific tissues and cell types.

Notice

- The *Pigs* gene is located on the Chr11. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of loxp insertion on gene transcription, RNA splicing and protein translation cannot be predicted at existing technological level.

Gene information (NCBI)

Pigs phosphatidylinositol glycan anchor biosynthesis, class S [Mus musculus (house mouse)]

Gene ID: 276846, updated on 17-Dec-2020

Summary



Official Symbol Pigs provided by [MGI](#)

Official Full Name phosphatidylinositol glycan anchor biosynthesis, class S provided by [MGI](#)

Primary source [MGI:MGI:2687325](#)

See related [Ensembl:ENSMUSG00000041958](#)

Gene type protein coding

RefSeq status PROVISIONAL

Organism [Mus musculus](#)

Lineage Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus

Also known as BC058979, Gm393, Gm689

Expression Ubiquitous expression in adrenal adult (RPKM 76.5), ovary adult (RPKM 51.5) and 28 other tissues [See more](#)


Orthologs [human](#) [all](#)

Transcript information (Ensembl)


The gene has 1 transcript, and the transcript is shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Pigs-201	ENSMUST00000048073.9	2501	555aa	Protein coding	CCDS25100		TSL:1 , GENCODE basic , APPRIS P1 ,


The strategy is based on the design of *Pigs-201* transcript, the transcription is shown below:



Genomic location distribution



Protein domain



If you have any questions, you are welcome to inquire.
Tel: 400-9660890

