

Cnksr1 Cas9-KO Strategy

Designer:

Huan Wang

Reviewer:

Huan Fan

Design Date:

2020-2-18

Project Overview

Project Name

Cnksr1

Project type

Cas9-KO

Strain background

C57BL/6JGpt

Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Cnksr1* gene. The schematic diagram is as follows:



- The *Cnksr1* gene has 5 transcripts. According to the structure of *Cnksr1* gene, exon2-exon14 of *Cnksr1-201* (ENSMUST00000030645.8) transcript is recommended as the knockout region. The region contains 1229bp coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Cnksr1* gene. The brief process is as follows: gRNA was transcribed in vitro. Cas9 and gRNA were microinjected into the fertilized eggs of C57BL/6JGpt mice. Fertilized eggs were transplanted to obtain positive F0 mice which were confirmed by PCR and sequencing. A stable F1 generation mouse model was obtained by mating positive F0 generation mice with C57BL/6JGpt mice.

- The *Cnksr1* gene is located on the Chr4. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This Strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of the gene knockout on gene transcription, RNA splicing and protein translation cannot be predicted at the existing technology level.

Gene information (NCBI)

Cnksr1 connector enhancer of kinase suppressor of Ras 1 [Mus musculus (house mouse)]

Gene ID: 194231, updated on 31-Jan-2019

Summary




Official Symbol	Cnksr1 provided by MGI
Official Full Name	connector enhancer of kinase suppressor of Ras 1 provided by MGI
Primary source	MGI:MGI:2670958
See related	Ensembl:ENSMUSG00000028841
Gene type	protein coding
RefSeq status	PROVISIONAL
Organism	Mus musculus
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	BC037758
Expression	Biased expression in large intestine adult (RPKM 14.7), colon adult (RPKM 13.4) and 14 other tissues See more
Orthologs	human all

Transcript information (Ensembl)


The gene has 5 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Cnksr1-201	ENSMUST00000030645.8	2468	700aa	Protein coding	CCDS38912	A2A9K7	TSL:1 GENCODE basic APPRIS P1
Cnksr1-203	ENSMUST00000144504.7	2395	No protein	Processed transcript	-	-	TSL:5
Cnksr1-202	ENSMUST00000124181.1	956	No protein	Processed transcript	-	-	TSL:5
Cnksr1-204	ENSMUST00000145998.7	450	No protein	Processed transcript	-	-	TSL:5
Cnksr1-205	ENSMUST00000146283.1	396	No protein	Processed transcript	-	-	TSL:3

The strategy is based on the design of *Cnksr1-201* transcript,The transcription is shown below




Genomic location distribution



Protein domain



集萃药康
GemPharmatech



If you have any questions, you are welcome to inquire.

Tel: 400-9660890

