

Cenph Cas9-KO Strategy

Designer:

Huan Wang

Design Date:

2019-7-24

Project Overview

Project Name

Cenph

Project type

Cas9-KO

Strain background

C57BL/6JGpt

Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Cenph* gene. The schematic diagram is as follows:



- The *Cenph* gene has 2 transcripts. According to the structure of *Cenph* gene, exon2 of *Cenph-201* (ENSMUST00000075550.3) transcript is recommended as the knockout region. The region contains 56bp coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Cenph* gene. The brief process is as follows: CRISPR/Cas9 system

- The *Cenph* gene is located on the Chr13. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This Strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of the gene knockout on gene transcription, RNA splicing and protein translation cannot be predicted at the existing technology level.

Gene information (NCBI)

Cenph centromere protein H [Mus musculus (house mouse)]

Gene ID: 26886, updated on 31-Jan-2019

Summary




Official Symbol	Cenph provided by MGI
Official Full Name	centromere protein H provided by MGI
Primary source	MGI:MGI:1349448
See related	Ensembl:ENSMUSG00000045273
Gene type	protein coding
RefSeq status	VALIDATED
Organism	Mus musculus
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	1700021I11Rik, 2410018A12Rik, 2610042E16Rik, 2810046K12Rik, AU044255, CENP-H, ENP
Expression	Biased expression in liver E14 (RPKM 15.8), CNS E11.5 (RPKM 12.8) and 9 other tissues See more
Orthologs	human all

Transcript information (Ensembl)


The gene has 2 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Cenph-202	ENSMUST00000136694.1	686	No protein	Retained intron	-	-	TSL:2
Cenph-201	ENSMUST00000075550.3	1251	241aa	Protein coding	CCDS26738	Q9QYM8	TSL:1 GENCODE basic APPRIS P1


The strategy is based on the design of *Cenph-201* transcript,The transcription is shown below




Genomic location distribution



Protein domain



Mouse phenotype description(MGI)



Phenotypes affected by the gene are marked in blue. Data quoted from MGI database(<http://www.informatics.jax.org/>).

If you have any questions, you are welcome to inquire.

Tel: 400-9660890

