

Erlec1 Cas9-KO Strategy

Designer:

Project Overview

Project Name

Erlec1

Project type

Cas9-KO

Strain background

C57BL/6JGpt

Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Erlec1* gene. The schematic diagram is as follows:



- The *Erlec1* gene has 5 transcripts. According to the structure of *Erlec1* gene, exon2-exon13 of *Erlec1-201* (ENSMUST00000073192.13) transcript is recommended as the knockout region. The region contains most of the coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Erlec1* gene. The brief process is as follows: CRISPR/Cas9 system

- The *Erlec1* gene is located on the Chr11. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This Strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of the gene knockout on gene transcription, RNA splicing and protein translation cannot be predicted at the existing technology level.

Gene information (NCBI)

Erlec1 endoplasmic reticulum lectin 1 [Mus musculus (house mouse)]

Gene ID: 66753, updated on 31-Jan-2019

Summary




Official Symbol	Erlec1 provided by MGI
Official Full Name	endoplasmic reticulum lectin 1 provided by MGI
Primary source	MGI:MGI:1914003
See related	Ensembl:ENSMUSG00000020311
Gene type	protein coding
RefSeq status	VALIDATED
Organism	Mus musculus
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	4933407N01Rik
Expression	Ubiquitous expression in cerebellum adult (RPKM 3.4), bladder adult (RPKM 3.3) and 28 other tissues See more
Orthologs	human all

Transcript information (Ensembl)


The gene has 5 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Erlec1-201	ENSMUST00000073192.13	2532	483aa	Protein coding	CCDS24509	Q8VEH8	TSL:1 GENCODE basic APPRIS P1
Erlec1-202	ENSMUST00000129593.8	977	238aa	Nonsense mediated decay	-	F6Z458	CDS 5' incomplete TSL:5
Erlec1-203	ENSMUST00000143126.2	366	80aa	Nonsense mediated decay	-	F7AQQ2	CDS 5' incomplete TSL:5
Erlec1-204	ENSMUST00000152770.1	532	No protein	Processed transcript	-	-	TSL:5
Erlec1-205	ENSMUST00000155304.1	703	No protein	Retained intron	-	-	TSL:3

The strategy is based on the design of *Erlec1-201* transcript,The transcription is shown below




Genomic location distribution



Protein domain



Mouse phenotype description(MGI)



Phenotypes affected by the gene are marked in blue. Data quoted from MGI database(<http://www.informatics.jax.org/>).

If you have any questions, you are welcome to inquire.

Tel: 400-9660890

