

Atp4b Cas9-CKO Strategy

Designer:

Huan Fan

Reviewer:

Huan Wang

Design Date:

2019-12-27

Project Overview

Project Name

Atp4b

Project type


Cas9-CKO

Strain background

C57BL/6JGpt

Conditional Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Atp4b* gene. The schematic diagram is as follows:



- The *Atp4b* gene has 2 transcripts. According to the structure of *Atp4b* gene, exon3-exon4 of *Atp4b-201* (ENSMUST00000033826.3) transcript is recommended as the knockout region. The region contains 314bp coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Atp4b* gene. The brief process is as follows: CRISPR/Cas9 system and Donor were microinjected into the fertilized eggs of C57BL/6JGpt mice. Fertilized eggs were transplanted to obtain positive F0 mice which were confirmed by PCR and sequencing. A stable F1 generation mouse model was obtained by mating positive F0 generation mice with C57BL/6JGpt mice.
- The flox mice will be knocked out after mating with mice expressing Cre recombinase, resulting in the loss of function of the target gene in specific tissues and cell types.

- According to the existing MGI data, Mice homozygous for disruptions in this gene are superficially normal but have hypertrophied stomach mucosa, neutral luminal pH in the stomach, and hypergastrinemia.
- The *Atp4b* gene is located on the Chr8. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This Strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of loxp insertion on gene transcription, RNA splicing and protein translation cannot be predicted at existing technological level.

Gene information (NCBI)

Atp4b ATPase, H⁺/K⁺ exchanging, beta polypeptide [Mus musculus (house mouse)]

Gene ID: 11945, updated on 31-Jan-2019

Summary




Official Symbol	Atp4b provided by MGI
Official Full Name	ATPase, H ⁺ /K ⁺ exchanging, beta polypeptide provided by MGI
Primary source	MGI:MGI:88114
See related	Ensembl:ENSMUSG00000031449
Gene type	protein coding
RefSeq status	VALIDATED
Organism	Mus musculus
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	AV080843
Expression	Restricted expression toward stomach adult (RPKM 1152.3) See more
Orthologs	human all

Transcript information (Ensembl)


The gene has 2 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Atp4b-201	ENSMUST00000033826.3	1378	294aa	Protein coding	CCDS22111	P50992 Q0VBB6	TSL:1 GENCODE basic APPRIS P1
Atp4b-202	ENSMUST00000210491.1	709	No protein	Retained intron	-	-	TSL:1

The strategy is based on the design of *Atp4b-201* transcript,The transcription is shown below




Genomic location distribution



Protein domain



Mouse phenotype description(MGI)



Phenotypes affected by the gene are marked in blue. Data quoted from MGI database(<http://www.informatics.jax.org/>).

According to the existing MGI data, Mice homozygous for disruptions in this gene are superficially normal but have hypertrophied stomach mucosa, neutral luminal pH in the stomach, and hypergastrinemia.

If you have any questions, you are welcome to inquire.

Tel: 400-9660890

