

Wdr34 Cas9-KO Strategy

Designer:

Huan Wang

Reviewer:

Huan Fan

Design Date:

2020-3-9

Project Overview

Project Name

Wdr34

Project type


Cas9-KO

Strain background

C57BL/6JGpt

Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Wdr34* gene. The schematic diagram is as follows:



- The *Wdr34* gene has 4 transcripts. According to the structure of *Wdr34* gene, exon2-exon3 of *Wdr34-201* (ENSMUST00000113711.2) transcript is recommended as the knockout region. The region contains 362bp coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Wdr34* gene. The brief process is as follows: CRISPR/Cas9 system

- The *Wdr34* gene is located on the Chr2. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This Strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of the gene knockout on gene transcription, RNA splicing and protein translation cannot be predicted at the existing technology level.

Gene information (NCBI)

Wdr34 WD repeat domain 34 [Mus musculus (house mouse)]

Gene ID: 71820, updated on 31-Jan-2019

Summary




Official Symbol	Wdr34 provided by MGI
Official Full Name	WD repeat domain 34 provided by MGI
Primary source	MGI:MGI:1919070
See related	Ensembl:ENSMUSG00000039715
Gene type	protein coding
RefSeq status	VALIDATED
Organism	Mus musculus
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	3200002I06Rik
Expression	Ubiquitous expression in ovary adult (RPKM 29.0), testis adult (RPKM 27.5) and 28 other tissues See more
Orthologs	human all

Transcript information (Ensembl)


The gene has 4 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Wdr34-201	ENSMUST00000113711.2	1814	537aa	Protein coding	CCDS15865	A2BE91	TSL:1 GENCODE basic APPRIS P1
Wdr34-202	ENSMUST00000139145.1	494	No protein	Retained intron	-	-	TSL:2
Wdr34-203	ENSMUST00000143892.1	454	No protein	Retained intron	-	-	TSL:2
Wdr34-204	ENSMUST00000146275.1	366	No protein	lncRNA	-	-	TSL:2


The strategy is based on the design of *Wdr34-201* transcript,The transcription is shown below



集萃药康
GemPharmatech



Protein domain



If you have any questions, you are welcome to inquire.

Tel: 400-9660890

