

Ucp3 Cas9-KO Strategy

Designer:

Huan Wang

Reviewer:

Huan Fan

Design Date:

2020-3-6

Project Overview

Project Name

Ucp3

Project type


Cas9-KO

Strain background

C57BL/6JGpt

Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Ucp3* gene. The schematic diagram is as follows:



- The *Ucp3* gene has 4 transcripts. According to the structure of *Ucp3* gene, exon2-exon6 of *Ucp3-201* (ENSMUST00000032958.13) transcript is recommended as the knockout region. The region contains start codon ATG. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Ucp3* gene. The brief process is as follows: CRISPR/Cas9 system

- According to the existing MGI data, Homozygous null mutants exhibit a lack of superoxide-induced uncoupling in skeletal muscle mitochondria, accompanied by increased reactive oxygen species formation.
- The *Ucp3* gene is located on the Chr7. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This Strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of the gene knockout on gene transcription, RNA splicing and protein translation cannot be predicted at the existing technology level.

Gene information (NCBI)

Ucp3 uncoupling protein 3 (mitochondrial, proton carrier) [Mus musculus (house mouse)]

Gene ID: 22229, updated on 9-Apr-2019

Summary




Official Symbol	Ucp3 provided by MGI
Official Full Name	uncoupling protein 3 (mitochondrial, proton carrier) provided by MGI
Primary source	MGI:MGI:1099787
See related	Ensembl:ENSMUSG00000032942
Gene type	protein coding
RefSeq status	VALIDATED
Organism	Mus musculus
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	AI645527, Slc25a9, UCP-3
Expression	Biased expression in mammary gland adult (RPKM 11.6), heart adult (RPKM 10.8) and 10 other tissues See more
Orthologs	human all

Transcript information (Ensembl)


The gene has 4 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Ucp3-201	ENSMUST00000032958.13	2448	308aa	Protein coding	CCDS21497	B2RTM2 P56501	TSL:1 GENCODE basic APPRIS P1
Ucp3-202	ENSMUST00000107059.1	1260	308aa	Protein coding	CCDS21497	B2RTM2 P56501	TSL:5 GENCODE basic APPRIS P1
Ucp3-203	ENSMUST00000128044.7	2129	No protein	Retained intron	-	-	TSL:1
Ucp3-204	ENSMUST00000133850.1	1615	No protein	Retained intron	-	-	TSL:1

The strategy is based on the design of *Ucp3-201* transcript,The transcription is shown below



Genomic location distribution



Protein domain

ENSMUSP00000032...

[Superfamily](#)

Mitochondrial carrier domain superfamily

[Prints](#)

Mitochondrial carrier UCP-like

[Pfam](#)

Mitochondrial substrate/solute carrier

[PROSITE profiles](#)

Mitochondrial substrate/solute carrier

[PANTHER](#)

PTHR45618:SF2

PTHR45618

[Gene3D](#)

Mitochondrial carrier domain superfamily

[All sequence SNPs/i...](#)

Sequence variants (dbSNP and all other sources)


[Variant Legend](#)

 synonymous variant

[Scale bar](#)

0 40 80 120 160 200 240 308

Mouse phenotype description(MGI)



Phenotypes affected by the gene are marked in blue. Data quoted from MGI database(<http://www.informatics.jax.org/>).

According to the existing MGI data, Homozygous null mutants exhibit a lack of superoxide-induced uncoupling in skeletal muscle mitochondria, accompanied by increased reactive oxygen species formation.

If you have any questions, you are welcome to inquire.

Tel: 400-9660890

