

Tnfrsf19 Cas9-KO Strategy

Designer:

Huan Fan

Design Date:

2019-8-23

Project Overview

Project Name

Tnfrsf19

Project type


Cas9-KO

Strain background

C57BL/6JGpt

Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Tnfrsf19* gene. The schematic diagram is as follows:



- The *Tnfrsf19* gene has 5 transcripts. According to the structure of *Tnfrsf19* gene, exon4 of *Tnfrsf19-201* (ENSMUST00000111234.9) transcript is recommended as the knockout region. The region contains 179bp coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Tnfrsf19* gene. The brief process is as follows: CRISPR/Cas9 system

- According to the existing MGI data, Homozygous mutant mice exhibit no obvious physical abnormalities or alterations in behavior, locomotion, or fecundity, however neurons are more resistant to the suppressive action of myelin inhibitors.
- The *Tnfrsf19* gene is located on the Chr14. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This Strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of the gene knockout on gene transcription, RNA splicing and protein translation cannot be predicted at the existing technology level.

Gene information (NCBI)

Tnfrsf19 tumor necrosis factor receptor superfamily, member 19 [Mus musculus (house mouse)]

Gene ID: 29820, updated on 3-Feb-2019

Summary




Official Symbol	Tnfrsf19 provided by MGI
Official Full Name	tumor necrosis factor receptor superfamily, member 19 provided by MGI
Primary source	MGI:MGI:1352474
See related	Ensembl:ENSMUSG00000060548
Gene type	protein coding
RefSeq status	VALIDATED
Organism	Mus musculus
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	AL023044, AW123854, TAJ, TAJ-ALPHA, TRADE, Troy
Expression	Broad expression in adrenal adult (RPKM 14.8), CNS E11.5 (RPKM 12.6) and 21 other tissues See more
Orthologs	human all

Transcript information (Ensembl)


The gene has 5 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Tnfrsf19-201	ENSMUST00000111234.9	4051	416aa	Protein coding	CCDS36941	Q9JLL3	TSL:1 GENCODE basic APPRIS P2
Tnfrsf19-202	ENSMUST00000111236.3	3958	416aa	Protein coding	CCDS36941	Q9JLL3	TSL:1 GENCODE basic APPRIS P2
Tnfrsf19-203	ENSMUST00000224371.1	1239	214aa	Protein coding	-	Q9JLL3	GENCODE basic APPRIS ALT2
Tnfrsf19-205	ENSMUST00000225730.1	744	214aa	Protein coding	-	Q9JLL3	GENCODE basic APPRIS ALT2
Tnfrsf19-204	ENSMUST00000225501.1	1087	No protein	Processed transcript	-	-	

The strategy is based on the design of *Tnfrsf19-201* transcript,The transcription is shown below




Genomic location distribution



Protein domain



Mouse phenotype description(MGI)



Phenotypes affected by the gene are marked in blue. Data quoted from MGI database(<http://www.informatics.jax.org/>).

According to the existing MGI data, Homozygous mutant mice exhibit no obvious physical abnormalities or alterations in behavior, locomotion, or fecundity, however neurons are more resistant to the suppressive action of myelin inhibitors.

If you have any questions, you are welcome to inquire.

Tel: 400-9660890

