

Irf6 Cas9-CKO Strategy

Designer:

Daohua Xu

Design Date:

2019-7-18

Project Overview

Project Name

Irf6

Project type


Cas9-CKO

Strain background

C57BL/6JGpt

Conditional Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Irf6* gene. The schematic diagram is as follows:



- The *Irf6* gene has 4 transcripts. According to the structure of *Irf6* gene, exon3-exon4 of *Irf6-201* (ENSMUST00000076521.6) transcript is recommended as the knockout region. The region contains 334bp coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Irf6* gene. The brief process is as follows: CRISPR/Cas9 system and Donor were microinjected into the fertilized eggs of C57BL/6JGpt mice. Fertilized eggs were transplanted to obtain positive F0 mice which were confirmed by PCR and sequencing. A stable F1 generation mouse model was obtained by mating positive F0 generation mice with C57BL/6JGpt mice.
- The flox mice will be knocked out after mating with mice expressing Cre recombinase, resulting in the loss of function of the target gene in specific tissues and cell types.

- According to the existing MGI data, Mice with mutations of *Irf6* display craniofacial, limb, and skin defects.
- The *Irf6* gene is located on the Chr1. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This Strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of loxp insertion on gene transcription, RNA splicing and protein translation cannot be predicted at existing technological level.

Gene information (NCBI)

Irf6 interferon regulatory factor 6 [Mus musculus (house mouse)]

Gene ID: 54139, updated on 23-Feb-2019

Summary




Official Symbol	Irf6 provided by MGI
Official Full Name	interferon regulatory factor 6 provided by MGI
Primary source	MGI:MGI:1859211
See related	Ensembl:ENSMUSG00000026638
Gene type	protein coding
RefSeq status	VALIDATED
Organism	Mus musculus
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	AI876454, E230028I05Rik, mirf6
Expression	Broad expression in large intestine adult (RPKM 19.7), colon adult (RPKM 15.5) and 16 other tissues See more
Orthologs	human all

Transcript information (Ensembl)


The gene has 4 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Irf6-201	ENSMUST00000076521.6	4113	467aa	Protein coding	CCDS15632	P97431	TSL:1 GENCODE basic APPRIS P1
Irf6-203	ENSMUST00000140066.2	2509	No protein	Retained intron	-	-	TSL:1
Irf6-202	ENSMUST00000138826.7	1760	No protein	Retained intron	-	-	TSL:1
Irf6-204	ENSMUST00000193772.1	341	No protein	Retained intron	-	-	TSL:3

The strategy is based on the design of *Irf6-201* transcript,The transcription is shown below




Genomic location distribution



Protein domain



Mouse phenotype description(MGI)



Phenotypes affected by the gene are marked in blue. Data quoted from MGI database(<http://www.informatics.jax.org/>).

According to the existing MGI data, Mice with mutations of *Irf6* display craniofacial, limb, and skin defects.

If you have any questions, you are welcome to inquire.

Tel: 400-9660890

