

Med19 Cas9-KO Strategy

Designer: Xueting Zhang

Project Overview

Project Name

Med19

Project type


Cas9-KO

Strain background

C57BL/6JGpt

Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Med19* gene. The schematic diagram is as follows:



- The *Med19* gene has 2 transcripts. According to the structure of *Med19* gene, exon2-exon4 of *Med19-201* (ENSMUST00000102645.3) transcript is recommended as the knockout region. The region contains 449bp coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Med19* gene. The brief process is as follows: gRNA was transcribed in vitro. Cas9 and gRNA were microinjected into the fertilized eggs of C57BL/6JGpt mice. Fertilized eggs were transplanted to obtain positive F0 mice which were confirmed by PCR and sequencing. A stable F1 generation mouse model was obtained by mating positive F0 generation mice with C57BL/6JGpt mice.

- The *Med19* gene is located on the Chr2. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This Strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of the gene knockout on gene transcription, RNA splicing and protein translation cannot be predicted at the existing technology level.

Gene information (NCBI)

Med19 mediator complex subunit 19 [*Mus musculus* (house mouse)]

Gene ID: 381379, updated on 31-Jan-2019

Summary

Official Symbol Med19 provided by [MGI](#)
Official Full Name mediator complex subunit 19 provided by [MGI](#)
Primary source [MGI:MGI:1914234](#)
See related [Ensembl:ENSMUSG00000027080](#)
Gene type protein coding
RefSeq status VALIDATED
Organism [Mus musculus](#)
Lineage Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as LCMR1; 2410018M14Rik; 3110040A13Rik
Expression Ubiquitous expression in CNS E14 (RPKM 11.2), whole brain E14.5 (RPKM 10.7) and 28 other tissues [See more](#)
Orthologs [human](#) [all](#)

Genomic context

Location: 2; 2 D

See Med19 in [Genome Data Viewer](#)

Exon count: 5


Annotation release	Status	Assembly	Chr	Location
106	current	GRCm38.p4 (GCF_000001635.24)	2	NC_000068.7 (84678402..84688215)
Build 37.2	previous assembly	MGSCv37 (GCF_000001635.18)	2	NC_000068.6 (84518559..84528372)

Transcript information (Ensembl)


The gene has 2 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Med19-201	ENSMUST00000102645.3	2059	244aa	Protein coding	CCDS16189	Q8C1S0	TSL:1 GENCODE basic APPRIS P1
Med19-202	ENSMUST00000141595.1	778	No protein	Processed transcript	-	-	TSL:3

The strategy is based on the design of *Med19-201* transcript,The transcription is shown below




Genomic location distribution



Protein domain



Mouse phenotype description(MGI)



Phenotypes affected by the gene are marked in blue. Data quoted from MGI database(<http://www.informatics.jax.org/>).

If you have any questions, you are welcome to inquire.

Tel: 400-9660890

