

Yeats4 Cas9-CKO Strategy

Designer: Yanhua Shen

Project Overview

Project Name

Yeats4

Project type


Cas9-CKO

Strain background

C57BL/6JGpt

Conditional Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Yeats4* gene. The schematic diagram is as follows:



- The *Yeats4* gene has 3 transcripts. According to the structure of *Yeats4* gene, exon3-exon6 of *Yeats4-201* (ENSMUST00000020382.7) transcript is recommended as the knockout region. The region contains 343bp coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Yeats4* gene. The brief process is as follows: CRISPR/Cas9 system and Donor were microinjected into the fertilized eggs of C57BL/6JGpt mice. Fertilized eggs were transplanted to obtain positive F0 mice which were confirmed by PCR and sequencing. A stable F1 generation mouse model was obtained by mating positive F0 generation mice with C57BL/6JGpt mice.
- The flox mice will be knocked out after mating with mice expressing Cre recombinase, resulting in the loss of function of the target gene in specific tissues and cell types.

Notice

- The *Yeats4* gene is located on the Chr10. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This Strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of loxp insertion on gene transcription, RNA splicing and protein translation cannot be predicted at existing technological level.

Gene information (NCBI)

Yeats4 YEATS domain containing 4 [Mus musculus (house mouse)]

Gene ID: 64050, updated on 31-Jan-2019

Summary




Official Symbol	Yeats4 provided by MGI
Official Full Name	YEATS domain containing 4 provided by MGI
Primary source	MGI:MGI:1927224
See related	Ensembl:ENSMUSG00000020171
Gene type	protein coding
RefSeq status	VALIDATED
Organism	Mus musculus
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	4930573H17Rik, B230215M10Rik, GAS41, NuBI-1
Expression	Ubiquitous expression in CNS E11.5 (RPKM 52.0), whole brain E14.5 (RPKM 41.1) and 28 other tissues See more
Orthologs	human all

Transcript information (Ensembl)


The gene has 3 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Yeats4-201	ENSMUST00000020382.7	1327	227aa	Protein coding	CCDS24191	Q9CR11	TSL:1 GENCODE basic APPRIS P1
Yeats4-202	ENSMUST00000218887.1	1339	95aa	Nonsense mediated decay	-	A0A1W2P7H0	TSL:1
Yeats4-203	ENSMUST00000219609.1	556	No protein	Retained intron	-	-	TSL:1


The strategy is based on the design of *Yeats4-201* transcript,The transcription is shown below



Genomic location distribution



Protein domain



If you have any questions, you are welcome to inquire.

Tel: 400-9660890

