

Fgf15 Cas9-KO Strategy

Designer:

Ruirui Zhang

Reviewer:

Huimin Su

Design Date:

2019-9-21

Project Overview

Project Name

Fgf15

Project type


Cas9-KO

Strain background

C57BL/6JGpt

Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Fgf15* gene. The schematic diagram is as follows:



- The *Fgf15* gene has 5 transcripts. According to the structure of *Fgf15* gene, exon1-exon2 of *Fgf15-201* (ENSMUST00000033389.6) transcript is recommended as the knockout region. The region contains start codon ATG. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Fgf15* gene. The brief process is as follows: CRISPR/Cas9 system

- According to the existing MGI data, Targeted inactivation of this gene leads to a severe underrepresentation of homozygotes at weaning as well as highly penetrant ventricular septal defects and malalignment of the aorta and pulmonary trunk. No abnormalities in otic induction or otic vesicle formation are observed.
- The KO region contains functional region of the *Gm26793* gene. Knockout the region may affect the function of *Gm26793* gene.
- The *Fgf15* gene is located on the Chr7. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This Strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of the gene knockout on gene transcription, RNA splicing and protein translation cannot be predicted at the existing technology level.

Gene information (NCBI)

Fgf15 fibroblast growth factor 15 [*Mus musculus* (house mouse)]

Gene ID: 14170, updated on 11-Sep-2019

Summary

Official Symbol	Fgf15 provided by MGI
Official Full Name	fibroblast growth factor 15 provided by MGI
Primary source	MGI:MGI:1096383
See related	Ensembl:ENSMUSG000000031073
Gene type	protein coding
RefSeq status	VALIDATED
Organism	Mus musculus
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	FGF19
Expression	Biased expression in CNS E11.5 (RPKM 7.5), large intestine adult (RPKM 6.2) and 3 other tissues See more
Orthologs	human all

Transcript information (Ensembl)


The gene has 5 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Fgf15-201	ENSMUST00000033389.6	1810	218aa	Protein coding	CCDS22053	Q35622 Q790L8	TSL:1 GENCODE basic APPRIS P1
Fgf15-204	ENSMUST00000207229.1	600	161aa	Protein coding	-	A0A140LIL5	CDS 3' incomplete TSL:3
Fgf15-203	ENSMUST00000207040.1	1753	No protein	Retained intron	-	-	TSL:NA
Fgf15-205	ENSMUST00000208515.1	809	No protein	Retained intron	-	-	TSL:2
Fgf15-202	ENSMUST00000207033.1	600	No protein	lncRNA	-	-	TSL:2


The strategy is based on the design of *Fgf15-201* transcript,The transcription is shown below




Genomic location distribution



Protein domain



Mouse phenotype description(MGI)



Phenotypes affected by the gene are marked in blue. Data quoted from MGI database(<http://www.informatics.jax.org/>).

According to the existing MGI data, Targeted inactivation of this gene leads to a severe underrepresentation of homozygotes at weaning as well as highly penetrant ventricular septal defects and malalignment of the aorta and pulmonary trunk. No abnormalities in otic induction or otic vesicle formation are observed.

If you have any questions, you are welcome to inquire.

Tel: 400-9660890

