

Mylip **Cas9-KO Strategy**

Designer:

Daohua Xu

Reviewer:

Huimin Su

Design Date:

2019-11-22

Project Overview

Project Name

Mylip

Project type

Cas9-KO

Strain background

C57BL/6JGpt

Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Myliip* gene. The schematic diagram is as follows:



- The *Myliip* gene has 2 transcripts. According to the structure of *Myliip* gene, exon3-exon6 of *Myliip-201* (ENSMUST00000038275.10) transcript is recommended as the knockout region. The region contains 970bp coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Myliip* gene. The brief process is as follows: CRISPR/Cas9 system

- The *Myliip* gene is located on the Chr13. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This Strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of the gene knockout on gene transcription, RNA splicing and protein translation cannot be predicted at the existing technology level.

Gene information (NCBI)

Mylip myosin regulatory light chain interacting protein [Mus musculus (house mouse)]

Gene ID: 218203, updated on 31-Jan-2019

Summary



Official Symbol	Mylip provided by MGI
Official Full Name	myosin regulatory light chain interacting protein provided by MGI
Primary source	MGI:MGI:2388271
See related	Ensembl:ENSMUSG00000038175
Gene type	protein coding
RefSeq status	PROVISIONAL
Organism	Mus musculus
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	9430057C20Rik, Mir
Expression	Ubiquitous expression in bladder adult (RPKM 23.3), mammary gland adult (RPKM 20.8) and 27 other tissues See more
Orthologs	human all

Transcript information (Ensembl)


The gene has 2 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Mylip-201	ENSMUST00000038275.10	3065	445aa	Protein coding	CCDS36648	Q8BM54	TSL:1 GENCODE basic APPRIS P1
Mylip-202	ENSMUST00000222178.1	1954	380aa	Protein coding	-	Q8BM54	TSL:1 GENCODE basic

The strategy is based on the design of *Mylip-201* transcript,The transcription is shown below



Genomic location distribution




Protein domain



集萃药康
GemPharmatech



Mouse phenotype description(MGI)



Phenotypes affected by the gene are marked in blue. Data quoted from MGI database(<http://www.informatics.jax.org/>).

If you have any questions, you are welcome to inquire.

Tel: 400-9660890

