

Zswim4 Cas9-KO Strategy

Designer:

Project Overview

Project Name

Zswim4

Project type


Cas9-KO

Strain background

C57BL/6JGpt

Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Zswim4* gene. The schematic diagram is as follows:



- The *Zswim4* gene has 1 transcript. According to the structure of *Zswim4* gene, exon2-exon7 of *Zswim4-201* (ENSMUST00000039480.6) transcript is recommended as the knockout region. The region contains 1378bp coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Zswim4* gene. The brief process is as follows: CRISPR/Cas9 system

- The *Zswim4* gene is located on the Chr8. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This Strategy is designed based on genetic information in existing databases. Due to the complexity of gene transcription and translation processes, all risks cannot be predicted under existing information.

Gene information (NCBI)

Zswim4 zinc finger SWIM-type containing 4 [Mus musculus (house mouse)]

Gene ID: 212168, updated on 31-Jan-2019

Summary




Official Symbol	Zswim4 provided by MGI
Official Full Name	zinc finger SWIM-type containing 4 provided by MGI
Primary source	MGI:MGI:2443726
See related	Ensembl:ENSMUSG00000035671
Gene type	protein coding
RefSeq status	PROVISIONAL
Organism	Mus musculus
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	D330002I24, E130119J17Rik
Expression	Ubiquitous expression in adrenal adult (RPKM 36.0), whole brain E14.5 (RPKM 21.2) and 26 other tissues See more
Orthologs	human all

Transcript information (Ensembl)


The gene has 1 transcript, and the transcript is shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Zswim4-201	ENSMUST00000039480.6	4625	1101aa	Protein coding	CCDS22471	Q8C7B8	TSL:1 GENCODE basic APPRIS P1


The strategy is based on the design of Zswim4-201 transcript, The transcription is shown below



Genomic location distribution



Protein domain



If you have any questions, you are welcome to inquire.

Tel: 400-9660890

