

Col9a2 Cas9-CKO Strategy

Designer:

JiaYu

Reviewer:

Xiaojing Li

Design Date:

2020-1-23

Project Overview

Project Name

Col9a2

Project type


Cas9-CKO

Strain background

C57BL/6JGpt

Conditional Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Col9a2* gene. The schematic diagram is as follows:



- The *Col9a2* gene has 3 transcripts. According to the structure of *Col9a2* gene, exon1-exon32 of *Col9a2-201* (ENSMUST00000030372.5) transcript is recommended as the knockout region. The region contains all of the coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Col9a2* gene. The brief process is as follows: CRISPR/Cas9 system and Donor were microinjected into the fertilized eggs of C57BL/6JGpt mice. Fertilized eggs were transplanted to obtain positive F0 mice which were confirmed by PCR and sequencing. A stable F1 generation mouse model was obtained by mating positive F0 generation mice with C57BL/6JGpt mice.
- The flox mice will be knocked out after mating with mice expressing Cre recombinase, resulting in the loss of function of the target gene in specific tissues and cell types.

Notice

- The *Col9a2* gene is located on the Chr4. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This Strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of loxp insertion on gene transcription, RNA splicing and protein translation cannot be predicted at existing technological level.

Gene information (NCBI)

Col9a2 collagen, type IX, alpha 2 [Mus musculus (house mouse)]

Gene ID: 12840, updated on 3-Feb-2019

Summary




Official Symbol	Col9a2 provided by MGI
Official Full Name	collagen, type IX, alpha 2 provided by MGI
Primary source	MGI:MGI:88466
See related	Ensembl:ENSMUSG00000028626
Gene type	protein coding
RefSeq status	VALIDATED
Organism	Mus musculus
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	AI427499, Col9a-2
Expression	Biased expression in limb E14.5 (RPKM 84.1), CNS E14 (RPKM 17.5) and 4 other tissues See more
Orthologs	human all

Transcript information (Ensembl)


The gene has 3 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Col9a2-201	ENSMUST00000030372.5	2924	688aa	Protein coding	CCDS18599	I7HJR1	TSL:1 GENCODE basic APPRIS P1
Col9a2-203	ENSMUST00000151987.1	492	No protein	Retained intron	-	-	TSL:3
Col9a2-202	ENSMUST00000140119.1	752	No protein	lncRNA	-	-	TSL:2


The strategy is based on the design of *Col9a2-201* transcript,The transcription is shown below




Genomic location distribution



Protein domain



Mouse phenotype description(MGI)



Phenotypes affected by the gene are marked in blue. Data quoted from MGI database(<http://www.informatics.jax.org/>).

If you have any questions, you are welcome to inquire.

Tel: 400-9660890

