

$Snx3$ Cas9-CKO Strategy

Designer: Xueting Zhang

Reviewer: Yanhua Shen

Date: 2020-02-19

Project Overview

Project Name

$Snx3$

Project type


Cas9-CKO

Strain background

C57BL/6JGpt

Conditional Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Snx3* gene. The schematic diagram is as follows:



Technical routes

- The *Snx3* gene has 3 transcripts. According to the structure of *Snx3* gene, exon2-exon3 of *Snx3-201* (ENSMUST00000019939.11) transcript is recommended as the knockout region. The region contains 221bp coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Snx3* gene. The brief process is as follows: CRISPR/Cas9 system and Donor were microinjected into the fertilized eggs of C57BL/6JGpt mice. Fertilized eggs were transplanted to obtain positive F0 mice which were confirmed by PCR and sequencing. A stable F1 generation mouse model was obtained by mating positive F0 generation mice with C57BL/6JGpt mice.
- The flox mice will be knocked out after mating with mice expressing Cre recombinase, resulting in the loss of function of the target gene in specific tissues and cell types.

Notice

- The *Snx3* gene is located on the Chr10. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This Strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of loxp insertion on gene transcription, RNA splicing and protein translation cannot be predicted at existing technological level.

Gene information (NCBI)

Snx3 sorting nexin 3 [*Mus musculus* (house mouse)]

(Gene ID: 54198, updated on 12-Nov-2019)

Summary

- Official Symbol** Snx3 provided by [MGI](#)
- Official Full Name** sorting nexin 3 provided by [MGI](#)
- Primary source** [MGI:MGI:1860188](#)
- See related** [Ensembl:ENSMUSG00000019804](#)
- Gene type** protein coding
- RefSeq status** VALIDATED
- Organism** [Mus musculus](#)
- Lineage** Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
- Also known as** SDP3
- Expression** Ubiquitous expression in placenta adult (RPKM 150.6), testis adult (RPKM 118.9) and 28 other tissues [See more](#)
- Orthologs** [human](#) [all](#)

Genomic context

Location: 10 B2; 10 22.89 cM [See Snx3 in Genome Data Viewer](#)

Exon count: 4


Annotation release	Status	Assembly	Chr	Location
108	current	GRCm38.p6 (GCF_000001635.26)	10	NC_000076.6 (42502054..42535369)
Build 37.2	previous assembly	MGSCv37 (GCF_000001635.18)	10	NC_000076.5 (42221860..42255175)

Transcript information (Ensembl)


The gene has 3 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Snx3-201	ENSMUST00000019939.11	1401	162aa	Protein coding	CCDS23812	Q78ZM0	TSL:1 GENCODE basic APPRIS P1
Snx3-202	ENSMUST00000105499.1	1057	130aa	Protein coding	-	D3Z789	TSL:5 GENCODE basic
Snx3-203	ENSMUST00000105500.7	753	140aa	Protein coding	-	D3Z6Z0	TSL:5 GENCODE basic


The strategy is based on the design of *Snx3-201* transcript,The transcription is shown below




Genomic location distribution



Protein domain



Mouse phenotype description(MGI)



Phenotypes affected by the gene are marked in blue. Data quoted from MGI database(<http://www.informatics.jax.org/>).

If you have any questions, you are welcome to inquire.

Tel: 400-9660890

