

Cdh18 Cas9-CKO Strategy

Designer:

Daohua Xu

Reviewer:

Huimin Su

Design Date:

2020-3-30

Project Overview

Project Name

Cdh18

Project type

Cas9-CKO

Strain background

C57BL/6JGpt

Conditional Knockout strategy


This model will use CRISPR/Cas9 technology to edit the *Cdh18* gene. The schematic diagram is as follows:

Donor and CRISPR/Cas9 System

Wild-type allele

Conditional KO allele

KO allele



- The *Cdh18* gene has 7 transcripts. According to the structure of *Cdh18* gene, exon3 of *Cdh18-203* (ENSMUST00000164787.7) transcript is recommended as the knockout region. The region contains 295bp coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Cdh18* gene. The brief process is as follows: CRISPR/Cas9 system and Donor were microinjected into the fertilized eggs of C57BL/6JGpt mice. Fertilized eggs were transplanted to obtain positive F0 mice which were confirmed by PCR and sequencing. A stable F1 generation mouse model was obtained by mating positive F0 generation mice with C57BL/6JGpt mice.
- The flox mice will be knocked out after mating with mice expressing Cre recombinase, resulting in the loss of function of the target gene in specific tissues and cell types.

Notice

- The *Cdh18* gene is located on the Chr15. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This Strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of loxp insertion on gene transcription, RNA splicing and protein translation cannot be predicted at existing technological level.

Gene information (NCBI)

Cdh18 cadherin 18 [*Mus musculus* (house mouse)]

Gene ID: 320865, updated on 13-Mar-2020

Summary


Official Symbol	Cdh18 provided by MGI
Official Full Name	cadherin 18 provided by MGI
Primary source	MGI:MGI:1344366
See related	Ensembl:ENSMUSG00000040420
Gene type	protein coding
RefSeq status	VALIDATED
Organism	<i>Mus musculus</i>
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	Cdh14l; B230220E17Rik; Tg(GFAP-APOE_i4)1Hol
Expression	Biased expression in cerebellum adult (RPKM 8.0), cortex adult (RPKM 2.6) and 5 other tissues See more
Orthologs	human all

Transcript information (Ensembl)


The gene has 7 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Cdh18-203	ENSMUST00000164787.7	2744	790aa	Protein coding	CCDS49585	E9Q9Q6	TSL:5 GENCODE basic APPRIS P1
Cdh18-207	ENSMUST00000226693.1	2165	575aa	Protein coding	-	A0A2I3BPS8	GENCODE basic
Cdh18-204	ENSMUST00000165614.8	2162	574aa	Protein coding	-	E9Q941	TSL:5 GENCODE basic
Cdh18-202	ENSMUST00000163361.7	1534	337aa	Protein coding	-	Q8BLB5	CDS 3' incomplete TSL:5
Cdh18-205	ENSMUST00000167623.1	1476	441aa	Protein coding	-	F6RSC2	CDS 5' incomplete TSL:5
Cdh18-201	ENSMUST00000040427.13	2320	No protein	Retained intron	-	-	TSL:1
Cdh18-206	ENSMUST00000170626.1	1819	No protein	Retained intron	-	-	TSL:1


The strategy is based on the design of *Cdh18-203* transcript,The transcription is shown below



Genomic location distribution



Protein domain



If you have any questions, you are welcome to inquire.

Tel: 400-9660890

