

Eif3f Cas9-CKO Strategy

Designer:

JiaYu

Reviewer:

Xiaojing Li

Design Date:

2020-2-24

Project Overview

Project Name

Eif3f

Project type


Cas9-CKO

Strain background

C57BL/6JGpt

Conditional Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Eif3f* gene. The schematic diagram is as follows:



Technical routes

- The *Eif3f* gene has 3 transcripts. According to the structure of *Eif3f* gene, exon2-exon4 of *Eif3f-201* (ENSMUST00000033342.6) transcript is recommended as the knockout region. The region contains 289bp coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Eif3f* gene. The brief process is as follows: CRISPR/Cas9 system and Donor were microinjected into the fertilized eggs of C57BL/6JGpt mice. Fertilized eggs were transplanted to obtain positive F0 mice which were confirmed by PCR and sequencing. A stable F1 generation mouse model was obtained by mating positive F0 generation mice with C57BL/6JGpt mice.
- The flox mice will be knocked out after mating with mice expressing Cre recombinase, resulting in the loss of function of the target gene in specific tissues and cell types.

Notice

- The *Eif3f* gene is located on the Chr7. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This Strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of loxp insertion on gene transcription, RNA splicing and protein translation cannot be predicted at existing technological level.

Gene information (NCBI)

Eif3f eukaryotic translation initiation factor 3, subunit F [Mus musculus (house mouse)]

Gene ID: 66085, updated on 31-Jan-2019

Summary




Official Symbol	Eif3f provided by MGI
Official Full Name	eukaryotic translation initiation factor 3, subunit F provided by MGI
Primary source	MGI:MGI:1913335
See related	Ensembl:ENSMUSG00000031029
Gene type	protein coding
RefSeq status	VALIDATED
Organism	Mus musculus
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	0610037M02Rik, AA409853, Eif3s5
Expression	Broad expression in adrenal adult (RPKM 509.7), stomach adult (RPKM 285.6) and 28 other tissues See more
Orthologs	human all

Transcript information (Ensembl)


The gene has 3 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Eif3f-201	ENSMUST00000033342.6	2514	361aa	Protein coding	CCDS21731	Q9DCH4	TSL:1 GENCODE basic APPRIS P1
Eif3f-202	ENSMUST00000147369.1	588	No protein	Retained intron	-	-	TSL:3
Eif3f-203	ENSMUST00000153596.1	439	No protein	Retained intron	-	-	TSL:1


The strategy is based on the design of *Eif3f-201* transcript,The transcription is shown below



Genomic location distribution



Protein domain



If you have any questions, you are welcome to inquire.

Tel: 400-9660890

