

Paqr5 Cas9-CKO Strategy

Designer:

Huimin Su

Reviewer:

Ruirui Zhang

Design Date:

2020-4-14

Project Overview

Project Name

Paqr5

Project type


Cas9-CKO

Strain background

C57BL/6JGpt

Conditional Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Paqr5* gene. The schematic diagram is as follows:



- The *Paqr5* gene has 3 transcripts. According to the structure of *Paqr5* gene, exon4 of *Paqr5-201* (ENSMUST00000034817.10) transcript is recommended as the knockout region. The region contains 206bp coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Paqr5* gene. The brief process is as follows: CRISPR/Cas9 system and Donor were microinjected into the fertilized eggs of C57BL/6JGpt mice. Fertilized eggs were transplanted to obtain positive F0 mice which were confirmed by PCR and sequencing. A stable F1 generation mouse model was obtained by mating positive F0 generation mice with C57BL/6JGpt mice.
- The flox mice will be knocked out after mating with mice expressing Cre recombinase, resulting in the loss of function of the target gene in specific tissues and cell types.

- The *Paqr5* gene is located on the Chr9. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of loxp insertion on gene transcription, RNA splicing and protein translation cannot be predicted at existing technological level.

Gene information (NCBI)

Paqr5 progesterin and adipoQ receptor family member V [Mus musculus (house mouse)]

Gene ID: 74090, updated on 20-Mar-2020

Summary




Official Symbol	Paqr5 provided by MGI
Official Full Name	progesterin and adipoQ receptor family member V provided by MGI
Primary source	MGI:MGI:1921340
See related	Ensembl:ENSMUSG00000032278
Gene type	protein coding
RefSeq status	PROVISIONAL
Organism	Mus musculus
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	0610010I15Rik, AV002411, BB115488, mPRg
Expression	Biased expression in testis adult (RPKM 97.4), kidney adult (RPKM 47.4) and 3 other tissues See more
Orthologs	human all

Transcript information (Ensembl)


The gene has 3 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Paqr5-201	ENSMUST00000034817.10	3566	330aa	Protein coding	CCDS23261	Q9DCU0	TSL:1 GENCODE basic APPRIS P2
Paqr5-202	ENSMUST00000113990.1	1073	316aa	Protein coding	-	Q497J6	TSL:1 GENCODE basic APPRIS ALT1
Paqr5-203	ENSMUST00000135050.1	3307	No protein	Processed transcript	-	-	TSL:1


The strategy is based on the design of *Paqr5-201* transcript,the transcription is shown below



Genomic location distribution



Protein domain



If you have any questions, you are welcome to inquire.

Tel: 400-9660890

