

Slc38a8 Cas9-CKO Strategy

Designer: Huimin Su

Reviewer: Ruiuri Zhang

Design Date: 2020-4-28

Project Overview

Project Name

Slc38a8

Project type


Cas9-CKO

Strain background

C57BL/6JGpt

Conditional Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Slc38a8* gene. The schematic diagram is as follows:



- The *Slc38a8* gene has 5 transcripts. According to the structure of *Slc38a8* gene, exon4-exon6 of *Slc38a8-201* (ENSMUST00000036748.14) transcript is recommended as the knockout region. The region contains 293bp coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Slc38a8* gene. The brief process is as follows: CRISPR/Cas9 system and Donor were microinjected into the fertilized eggs of C57BL/6JGpt mice. Fertilized eggs were transplanted to obtain positive F0 mice which were confirmed by PCR and sequencing. A stable F1 generation mouse model was obtained by mating positive F0 generation mice with C57BL/6JGpt mice.
- The flox mice will be knocked out after mating with mice expressing Cre recombinase, resulting in the loss of function of the target gene in specific tissues and cell types.

Notice

- The *Slc38a8* gene is located on the Chr8. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of loxp insertion on gene transcription, RNA splicing and protein translation cannot be predicted at existing technological level.

Gene information (NCBI)

Slc38a8 solute carrier family 38, member 8 [*Mus musculus* (house mouse)]

Gene ID: 234788, updated on 13-Mar-2020

Summary




Official Symbol	Slc38a8 provided by MGI
Official Full Name	solute carrier family 38, member 8 provided by MGI
Primary source	MGI:MGI:2685433
See related	Ensembl:ENSMUSG00000034224
Gene type	protein coding
RefSeq status	VALIDATED
Organism	Mus musculus
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	Gm587
Expression	Low expression observed in reference dataset See more
Orthologs	human all

Transcript information (Ensembl)


The gene has 5 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Slc38a8-201	ENSMUST00000036748.14	2623	432aa	Protein coding	CCDS22706	Q5HZH7	TSL:1 GENCODE basic APPRIS P1
Slc38a8-205	ENSMUST00000138061.7	724	206aa	Protein coding	-	D3Z7F9	CDS 3' incomplete TSL:3
Slc38a8-204	ENSMUST00000133821.1	719	165aa	Protein coding	-	D3YUL2	CDS 3' incomplete TSL:3
Slc38a8-203	ENSMUST00000132838.7	869	No protein	Processed transcript	-	-	TSL:3
Slc38a8-202	ENSMUST00000125594.1	899	No protein	Retained intron	-	-	TSL:5


The strategy is based on the design of *Slc38a8-201* transcript, the transcription is shown below



Genomic location distribution



Protein domain



If you have any questions, you are welcome to inquire.

Tel: 400-9660890

