

Txnrd2 Cas9-CKO Strategy

Designer: Xueting Zhang

Reviewer: Yanhua Shen

Design Date: 2020-4-23

Project Overview



Project Name

Txnrd2

Project type

Cas9-CKO


Strain background

C57BL/6JGpt

Conditional Knockout strategy



This model will use CRISPR/Cas9 technology to edit the *Txnrd2* gene. The schematic diagram is as follows:



Technical routes



- The *Txnrd2* gene has 16 transcripts. According to the structure of *Txnrd2* gene, exon3-exon8 of *Txnrd2-215*(ENSMUST00000206151.1) transcript is recommended as the knockout region. The region contains 490bp coding sequence.

 Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Txnrd2* gene. The brief process is as follows:CRISPR/Cas9 system and Donor were microinjected into the fertilized eggs of C57BL/6JGpt mice. Fertilized eggs were transplanted to obtain positive F0 mice which were confirmed by PCR and sequencing. A stable F1 generation mouse model was obtained by mating positive F0 generation mice with C57BL/6JGpt mice.
- The flox mice will be knocked out after mating with mice expressing Cre recombinase, resulting in the loss of function of the target gene in specific tissues and cell types.

Notice



- ➤ According to the existing MGI data, mice homozygous for a knock-out allele die at e13 due to severe anemia and growth retardation, resulting from perturbed cardiac development and augmented apoptosis of hematopoietic cells.
- ➤ The effect on transcript *Txnrd2*-207 is unknown.
- ➤ Transcript *Txnrd2*-204&209&210&211 may not be affected.
- The *Txnrd2* gene is located on the Chr16. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of loxp insertion on gene transcription, RNA splicing and protein translation cannot be predicted at existing technological level.

Gene information (NCBI)



Txnrd2 thioredoxin reductase 2 [Mus musculus (house mouse)]

Gene ID: 26462, updated on 13-Mar-2020

Summary

△ ?

Official Symbol Txnrd2 provided by MGI

Official Full Name thioredoxin reductase 2 provided by MGI

Primary source MGI:MGI:1347023

See related Ensembl: ENSMUSG00000075704

Gene type protein coding
RefSeq status REVIEWED
Organism Mus musculus

Lineage Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha;

Muroidea; Muridae; Murinae; Mus; Mus

Also known as AA118373, ESTM573010, TGR, Tr3, Trxr2, Trxrd2

Summary The protein encoded by this gene belongs to the pyridine nucleotide-disulfide oxidoreductase family, and is a member of the thioredoxin

(Trx) system. Three thioredoxin reductase (TrxR) isozymes are found in mammals. TrxRs are selenocysteine-containing flavoenzymes, which reduce thioredoxins, as well as other substrates, and play a key role in redox homoeostasis. This gene encodes a mitochondrial form important for scavenging reactive oxygen species in mitochondria. It functions as a homodimer containing FAD, and selenocysteine (Sec) at the active site. Sec is encoded by UGA codon that normally signals translation termination. The 3' UTRs of selenoprotein mRNAs contain a conserved stem-loop structure, the Sec insertion sequence (SECIS) element, which is necessary for the recognition of UGA as a Sec codon rather than as a stop signal. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided

by RefSeq, Jun 2017]

Expression Ubiquitous expression in liver E14.5 (RPKM 41.1), liver E14 (RPKM 39.1) and 27 other tissuesSee more

Orthologs human all


Transcript information (Ensembl)



The gene has 16 transcripts, all transcripts are shown below:


Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
ENSMUST00000206151.1	2874	<u>527aa</u>	Protein coding	CCDS79434	A0A0M3HEQ0	TSL:5 GENCODE basic APPRIS is a system to annotate alternatively spliced transcripts based on a range of computational methods to identify the most functionally important transcript(s) of a gene. APPRIS P2
ENSMUST00000115606.7	1905	<u>527aa</u>	Protein coding	CCDS79434	A0A0M3HEQ0	TSL:1 GENCODE basic APPRIS is a system to annotate alternatively spliced transcripts based on a range of computational methods to identify the most functionally important transcript(s) of a gene. APPRIS P2
ENSMUST00000177856.7	1878	522aa	Protein coding	1920	J3QMQ8	TSL:5 GENCODE basic APPRIS is a system to annotate alternatively spliced transcripts based on a range of computational methods to identify the most functionally important transcript(s) of a gene. APPRIS ALT
ENSMUST00000205679.1	1866	505aa	Protein coding	1/27	A0A0U1RQA4	TSL:5 GENCODE basic
ENSMUST00000206606.1	1809	496aa	Protein coding		A0A0U1RPV8	TSL:1 GENCODE basic
ENSMUST00000178093.2	1785	491aa	Protein coding		J3QMN4	TSL:5 GENCODE basic
ENSMUST00000115604.7	1235	359aa	Protein coding	(2)	A0A0M3HEP9	TSL:1 GENCODE basic
ENSMUST00000138310.2	640	96aa	Protein coding	121	A0A0U1RPY5	CDS 5' incomplete TSL:3
ENSMUST00000144233.2	593	<u>165aa</u>	Protein coding	-	A0A0U1RPC6	CDS 3' incomplete TSL:5
ENSMUST00000131303.7	532	<u>123aa</u>	Protein coding		A0A0U1RPS1	CDS 5' incomplete TSL:2
ENSMUST00000126778.7	491	101aa	Protein coding	(40)	Q9D8I4	TSL:1 GENCODE basic
ENSMUST00000115605.7	1800	466aa	Nonsense mediated decay	121	D3Z0K8	TSL:5
ENSMUST00000149937.7	720	No protein	Processed transcript		-	TSL:5
ENSMUST00000153457.7	648	No protein	Processed transcript	(#)	-	TSL:5
ENSMUST00000156604.2	403	No protein	Processed transcript	(40)		TSL:5
ENSMUST00000126330.7	266	No protein	Processed transcript	100		TSL:5
	ENSMUST00000206151.1 ENSMUST00000115606.7 ENSMUST00000177856.7 ENSMUST00000205679.1 ENSMUST00000178093.2 ENSMUST00000115604.7 ENSMUST00000138310.2 ENSMUST00000144233.2 ENSMUST00000144233.2 ENSMUST000001478.7 ENSMUST00000115605.7 ENSMUST00000149937.7 ENSMUST00000153457.7 ENSMUST00000156604.2	ENSMUST00000206151.1 2874 ENSMUST00000115606.7 1905 ENSMUST00000177856.7 1878 ENSMUST00000205679.1 1866 ENSMUST00000178093.2 1785 ENSMUST00000115604.7 1235 ENSMUST00000115604.7 1235 ENSMUST000001138310.2 640 ENSMUST00000134233.2 593 ENSMUST00000134233.2 593 ENSMUST00000156778.7 491 ENSMUST00000115605.7 1800 ENSMUST00000115605.7 720 ENSMUST00000153457.7 648 ENSMUST00000156604.2 403	ENSMUST0000015606.7 1905 527aa ENSMUST00000115606.7 1905 527aa ENSMUST00000177856.7 1878 522aa ENSMUST00000205679.1 1866 505aa ENSMUST00000206606.1 1809 496aa ENSMUST00000115604.7 1235 359aa ENSMUST00000113604.7 1235 359aa ENSMUST00000138310.2 640 96aa ENSMUST0000013423.2 593 165aa ENSMUST0000013423.2 593 165aa ENSMUST00000144233.2 593 105aa ENSMUST0000015605.7 1800 466aa ENSMUST00000115605.7 1800 466aa ENSMUST00000153457.7 648 No protein ENSMUST00000156604.2 403 No protein	ENSMUST00000206151.1 2874 527aa Protein coding ENSMUST00000115606.7 1905 527aa Protein coding ENSMUST00000177856.7 1878 522aa Protein coding ENSMUST00000205679.1 1866 505aa Protein coding ENSMUST00000206606.1 1809 496aa Protein coding ENSMUST00000178093.2 1785 491aa Protein coding ENSMUST00000115604.7 1235 359aa Protein coding ENSMUST00000138310.2 640 96aa Protein coding ENSMUST00000144233.2 593 165aa Protein coding ENSMUST00000131303.7 532 123aa Protein coding ENSMUST00000126778.7 491 101aa Protein coding ENSMUST0000015605.7 1800 466aa Nonsense mediated decay ENSMUST00000153457.7 648 No protein Processed transcript ENSMUST00000156604.2 403 No protein Processed transcript	ENSMUST00000206151.1 2874 527aa Protein coding CCDS79434 ENSMUST00000115606.7 1905 527aa Protein coding CCDS79434 ENSMUST00000177856.7 1878 522aa Protein coding - ENSMUST00000205679.1 1866 505aa Protein coding - ENSMUST000001206606.1 1809 496aa Protein coding - ENSMUST00000118093.2 1785 491aa Protein coding - ENSMUST00000138310.2 640 96aa Protein coding - ENSMUST00000144233.2 593 165aa Protein coding - ENSMUST00000131303.7 532 123aa Protein coding - ENSMUST00000126778.7 491 101aa Protein coding - ENSMUST00000149937.7 720 No protein Processed transcript - ENSMUST00000153457.7 648 No protein Processed transcript - ENSMUST00000156604.2 403 No protein Processed transcript -	ENSMUST00000206151.1 2874 527aa Protein coding CCDS79434 A0A0M3HEQ0 ENSMUST00000115606.7 1905 527aa Protein coding CCDS79434 A0A0M3HEQ0 ENSMUST00000177856.7 1878 522aa Protein coding - J3QMQ8 ENSMUST00000205679.1 1866 505aa Protein coding - A0A0U1RQA4 ENSMUST00000178093.2 1785 491aa Protein coding - A0A0U1RPV8 ENSMUST00000115604.7 1235 359aa Protein coding - A0A0M3HEP9 ENSMUST00000138310.2 640 96aa Protein coding - A0A0U1RPV5 ENSMUST00000144233.2 593 165aa Protein coding - A0A0U1RPC6 ENSMUST00000131303.7 532 123aa Protein coding - A0A0U1RPS1 ENSMUST00000126778.7 491 101aa Protein coding - A0A0U1RPS1 ENSMUST0000015605.7 1800 466aa Nonsense mediated decay - D3Z0K8 ENSMUST00000153457.7<

The strategy is based on the design of *Txnrd2-215* transcript, the transcription is shown below:




Genomic location distribution






Protein domain





Mouse phenotype description(MGI)





Phenotypes affected by the gene are marked in blue.Data quoted from MGI database(http://www.informatics.jax.org/).

According to the existing MGI data, Mice homozygous for a knock-out allele die at E13 due to severe anemia and growth retardation, resulting from perturbed cardiac development and augmented apoptosis of hematopoietic cells.



If you have any questions, you are welcome to inquire. Tel: 400-9660890





