

Dnaaf5 Cas9-CKO Strategy

Designer: Xiaojing Li

Reviewer: Jia Yu

Design Date: 2020-4-20

Project Overview

Project Name

Dnaaf5

Project type


Cas9-CKO

Strain background

C57BL/6JGpt

Conditional Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Dnaaf5* gene. The schematic diagram is as follows:



- The *Dnaaf5* gene has 3 transcripts. According to the structure of *Dnaaf5* gene, exon5-exon6 of *Dnaaf5*-201 (ENSMUST00000026975.10) transcript is recommended as the knockout region. The region contains 446bp coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Dnaaf5* gene. The brief process is as follows: CRISPR/Cas9 system and Donor were microinjected into the fertilized eggs of C57BL/6JGpt mice. Fertilized eggs were transplanted to obtain positive F0 mice which were confirmed by PCR and sequencing. A stable F1 generation mouse model was obtained by mating positive F0 generation mice with C57BL/6JGpt mice.
- The flox mice will be knocked out after mating with mice expressing Cre recombinase, resulting in the loss of function of the target gene in specific tissues and cell types.

Notice

- The *Dnaaf5* gene is located on the Chr5. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This Strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of loxp insertion on gene transcription, RNA splicing and protein translation cannot be predicted at existing technological level.

Gene information (NCBI)

Dnaaf5 dynein, axonemal assembly factor 5 [Mus musculus (house mouse)]

Gene ID: 433956, updated on 13-Mar-2020

Summary




Official Symbol	Dnaaf5 provided by MGI
Official Full Name	dynein, axonemal assembly factor 5 provided by MGI
Primary source	MGI:MGI:3616079
See related	Ensembl:ENSMUSG00000025857
Gene type	protein coding
RefSeq status	VALIDATED
Organism	Mus musculus
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	BC053401, C76907, Heatr2
Expression	Ubiquitous expression in adrenal adult (RPKM 20.8), ovary adult (RPKM 20.2) and 28 other tissues See more
Orthologs	human all

Transcript information (Ensembl)


The gene has 3 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Dnaaf5-201	ENSMUST00000026975.10	3450	853aa	Protein coding	CCDS39349	B9EJR8	TSL:1 GENCODE basic APPRIS P1
Dnaaf5-203	ENSMUST00000196441.1	2426	558aa	Protein coding	-	B2RY24	TSL:1 GENCODE basic
Dnaaf5-202	ENSMUST00000127542.1	470	No protein	Processed transcript	-	-	TSL:3


The strategy is based on the design of *Dnaaf5-201* transcript,The transcription is shown below



Genomic location distribution



Protein domain



If you have any questions, you are welcome to inquire.

Tel: 400-9660890

