

# *Olfr24* Cas9-CKO Strategy

Designer: Yanhua Shen

Reviewer: Xueting Zhang

Design Date: 2020-2-28

# Project Overview

**Project Name**

***Olfr24***

**Project type**


**Cas9-CKO**

**Strain background**

**C57BL/6JGpt**

# Conditional Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Olfr24* gene. The schematic diagram is as follows:



# Technical routes

- The *Olfcr24* gene has 3 transcripts. According to the structure of *Olfcr24* gene, exon2 of *Olfcr24-203* (ENSMUST00000216754.2) transcript is recommended as the knockout region. The region contains all of the coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Olfcr24* gene. The brief process is as follows: CRISPR/Cas9 system and Donor were microinjected into the fertilized eggs of C57BL/6JGpt mice. Fertilized eggs were transplanted to obtain positive F0 mice which were confirmed by PCR and sequencing. A stable F1 generation mouse model was obtained by mating positive F0 generation mice with C57BL/6JGpt mice.
- The flox mice will be knocked out after mating with mice expressing Cre recombinase, resulting in the loss of function of the target gene in specific tissues and cell types.

# Notice

- The *Olfir24* gene is located on the Chr9. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This Strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of loxp insertion on gene transcription, RNA splicing and protein translation cannot be predicted at existing technological level.



# Gene information (NCBI)

## Olfr24 olfactory receptor 24 [ *Mus musculus* (house mouse) ]

Gene ID: 18322, updated on 10-Oct-2019

### Summary

Official Symbol	Olfr24 provided by MGI
Official Full Name	olfactory receptor 24 provided by MGI
Primary source	MGI:MGI:109311
See related	Ensembl:ENSMUSG00000054141
Gene type	protein coding
RefSeq status	PROVISIONAL
Organism	<i>Mus musculus</i>
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	MTPCR51; MOR132-1
Summary	<p>Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]</p>
Orthologs	human <a href="#">all</a>

### Genomic context

Location: 9; 9 A2

Exon count: 1


See Olfr24 in [Genome Data Viewer](#)

# Transcript information (Ensembl)


The gene has 3 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Olfr24-203	<a href="#">ENSMUST00000216754.2</a>	4327	<a href="#">313aa</a>	Protein coding	<a href="#">CCDS22846</a>	<a href="#">Q8VFM9</a>	TSL:2 GENCODE basic APPRIS P1
Olfr24-201	<a href="#">ENSMUST00000066997.3</a>	1157	<a href="#">313aa</a>	Protein coding	<a href="#">CCDS22846</a>	<a href="#">Q8VFM9</a>	TSL:NA GENCODE basic APPRIS P1
Olfr24-202	<a href="#">ENSMUST00000214366.1</a>	3505	No protein	lncRNA	-	-	TSL:5

The strategy is based on the design of *Olfr24-203* transcript,The transcription is shown below




# Genomic location distribution





# Protein domain



If you have any questions, you are welcome to inquire.

Tel: 400-9660890

