

***Kdelr2* Cas9-KO Strategy**

Designer:

Huimin Su

Reviewer:

Ruirui Zhang

Design Date:

2019/9/11

Project Overview

Project Name

Kdelr2

Project type

Cas9-KO

Strain background

C57BL/6JGpt

Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Kdelr2* gene. The schematic diagram is as follows:



- The *Kdelr2* gene has 1 transcript. According to the structure of *Kdelr2* gene, exon2 of *Kdelr2-201* (ENSMUST00000110731.3) transcript is recommended as the knockout region. The region contains 101bp coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Kdelr2* gene. The brief process is as follows: CRISPR/Cas9 system

- The *Kdelr2* gene is located on the Chr5. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This Strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of the gene knockout on gene transcription, RNA splicing and protein translation cannot be predicted at the existing technology level.

Gene information (NCBI)

Kdelr2 KDEL (Lys-Asp-Glu-Leu) endoplasmic reticulum protein retention receptor 2 [*Mus musculus* (house mouse)]

Gene ID: 66913, updated on 14-Aug-2019

Summary


Official Symbol	Kdelr2 provided by MGI
Official Full Name	KDEL (Lys-Asp-Glu-Leu) endoplasmic reticulum protein retention receptor 2 provided by MGI
Primary source	MGI:MGI:1914163
See related	Ensembl:ENSMUSG00000079111
Gene type	protein coding
RefSeq status	VALIDATED
Organism	Mus musculus
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	1110007A14Rik
Expression	Ubiquitous expression in duodenum adult (RPKM 116.2), colon adult (RPKM 108.7) and 28 other tissues See more
Orthologs	human all

Transcript information (Ensembl)


The gene has 1 transcript, and the transcript is shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Kdelr2-201	ENSMUST00000110731.3	1829	212aa	Protein coding	CCDS39369	Q9CQM2	TSL:1 Gencode basic APPRIS P1


The strategy is based on the design of *Kdelr2-201* transcript, The transcription is shown below



Genomic location distribution



Protein domain



If you have any questions, you are welcome to inquire.

Tel: 400-9660890

