

***Grik4* Cas9-KO Strategy**

Designer: Xueting Zhang

Reviewer: Yanhua Shen

Date: 2020-1-14

Project Overview

Project Name

Grik4

Project type

Cas9-KO

Strain background

C57BL/6JGpt

Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Grik4* gene. The schematic diagram is as follows:



- The *Grik4* gene has 2 transcripts. According to the structure of *Grik4* gene, exon4-exon5 of *Grik4-202* (ENSMUST00000114865.7) transcript is recommended as the knockout region. The region contains 263bp coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Grik4* gene. The brief process is as follows: CRISPR/Cas9 system

- According to the existing MGI data, Mice homozygous for a knock-out allele exhibit reduced GYKI-resistant excitatory postsynaptic current.
- The *Grik4* gene is located on the Chr9. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This Strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of the gene knockout on gene transcription, RNA splicing and protein translation cannot be predicted at the existing technology level.

Gene information (NCBI)

Grik4 glutamate receptor, ionotropic, kainate 4 [*Mus musculus* (house mouse)]

Gene ID: 110637, updated on 10-Oct-2019

Summary

Official Symbol	Grik4 provided by MGI
Official Full Name	glutamate receptor, ionotropic, kainate 4 provided by MGI
Primary source	MGI:MGI:95817
See related	Ensembl:ENSMUSG00000032017
Gene type	protein coding
RefSeq status	VALIDATED
Organism	<i>Mus musculus</i>
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	KA1; KA-1; GluK4; GluRgamma1; 6330551K01Rik
Expression	Broad expression in cortex adult (RPKM 1.7), frontal lobe adult (RPKM 1.6) and 18 other tissues See more
Orthologs	human all

Genomic context

Location: 9 A5.1; 9 23.8 cM

See Grik4 in [Genome Data Viewer](#)

Exon count: 27


Annotation release	Status	Assembly	Chr	Location
108	current	GRCm38.p6 (GCF_000001635.26)	9	NC_000075.6 (42518182..42945000, complement)
Build 37.2	previous assembly	MGSCv37 (GCF_000001635.18)	9	NC_000075.5 (42328519..42752454, complement)

Transcript information (Ensembl)


The gene has 2 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Grik4-202	ENSMUST00000114865.7	5887	956aa	Protein coding	CCDS23089	Q8BMF5	TSL:1 GENCODE basic APPRIS P1
Grik4-201	ENSMUST00000034515.6	3413	956aa	Protein coding	CCDS23089	Q8BMF5	TSL:5 GENCODE basic APPRIS P1


The strategy is based on the design of *Grik4-202* transcript,The transcription is shown below




Genomic location distribution



Protein domain



Mouse phenotype description(MGI)



Phenotypes affected by the gene are marked in blue. Data quoted from MGI database(<http://www.informatics.jax.org/>).

According to the existing MGI data, Mice homozygous for a knock-out allele exhibit reduced GYKI-resistant excitatory postsynaptic current.

If you have any questions, you are welcome to inquire.

Tel: 400-9660890

