

Matn3 Cas9-KO Strategy

Designer: Xueting Zhang

Reviewer: Yanhua Shen

Date: 2020-02-06

Project Overview

Project Name

Matn3

Project type

Cas9-KO

Strain background

C57BL/6JGpt

Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Matn3* gene. The schematic diagram is as follows:



- The *Matn3* gene has 1 transcript. According to the structure of *Matn3* gene, exon2-exon7 of *Matn3-201* (ENSMUST00000020899.4) transcript is recommended as the knockout region. The region contains most of the coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Matn3* gene. The brief process is as follows: CRISPR/Cas9 system

- According to the existing MGI data, Mice homozygous for disruptions in this gene display a normal phenotype.
- The knockout region is near to the N-terminal of *Wdr35* gene, this strategy may influence the regulatory function of the N-terminal of *Wdr35* gene.
- The *Matn3* gene is located on the Chr12. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This Strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of the gene knockout on gene transcription, RNA splicing and protein translation cannot be predicted at the existing technology level.

Gene information (NCBI)

Matn3 matrilin 3 [*Mus musculus* (house mouse)]

Gene ID: 17182, updated on 21-Aug-2019

Summary

- Official Symbol** Matn3 provided by [MGI](#)
- Official Full Name** matrilin 3 provided by [MGI](#)
- Primary source** [MGI:MGI:1328350](#)
- See related** [Ensembl:ENSMUSG00000020583](#)
- Gene type** protein coding
- RefSeq status** VALIDATED
- Organism** [Mus musculus](#)
- Lineage** Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
- Also known as** AV009181
- Expression** Restricted expression toward limb E14.5 (RPKM 28.9) [See more](#)
- Orthologs** [human](#) [all](#)

Genomic context

Location: 12 A1.1; 12 3.96 cM

See Matn3 in [Genome Data Viewer](#)

Exon count: 8

Annotation release	Status	Assembly	Chr	Location
108	current	GRCm38.p6 (GCF_000001635.26)	12	NC_000078.6 (8947929..8972028)
Build 37.2	previous assembly	MGSCv37 (GCF_000001635.18)	12	NC_000078.5 (8954735..8978299)

Transcript information (Ensembl)


The gene has 1 transcript, and the transcript is shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Matn3-201	ENSMUST00000020899.4	3014	481aa	Protein coding	CCDS25806	Q35701	TSL:1 GENCODE basic APPRIS P1


The strategy is based on the design of *Matn3-201* transcript, The transcription is shown below




Genomic location distribution



Protein domain



Mouse phenotype description(MGI)



Phenotypes affected by the gene are marked in blue. Data quoted from MGI database(<http://www.informatics.jax.org/>).

According to the existing MGI data, Mice homozygous for disruptions in this gene display a normal phenotype.

If you have any questions, you are welcome to inquire.

Tel: 400-9660890

