

Prrg4 Cas9-KO Strategy

Designer: Xueting Zhang

Reviewer: Daohua Xu

Design Date: 2020-5-14

Project Overview

Project Name

Prrg4

Project type


Cas9-KO

Strain background

C57BL/6JGpt

Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Prrg4* gene. The schematic diagram is as follows:



- The *Prrg4* gene has 2 transcripts. According to the structure of *Prrg4* gene, exon3 of *Prrg4-201* (ENSMUST00000028593.10) transcript is recommended as the knockout region. The region contains 164bp coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Prrg4* gene. The brief process is as follows: CRISPR/Cas9 system v

- The *Prrg4* gene is located on the Chr2. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of the gene knockout on gene transcription, RNA splicing and protein translation cannot be predicted at the existing technology level.

Gene information (NCBI)

Prrg4 proline rich Gla (G-carboxyglutamic acid) 4 (transmembrane) [Mus musculus (house mouse)]

Gene ID: 228413, updated on 13-Mar-2020

Summary




Official Symbol	Prrg4 provided by MGI
Official Full Name	proline rich Gla (G-carboxyglutamic acid) 4 (transmembrane) provided by MGI
Primary source	MGI:MGI:2442211
See related	Ensembl:ENSMUSG000000027171
Gene type	protein coding
RefSeq status	VALIDATED
Organism	Mus musculus
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	993011118Rik, TMG4
Expression	Biased expression in placenta adult (RPKM 2.9), bladder adult (RPKM 1.7) and 14 other tissues See more
Orthologs	human all

Transcript information (Ensembl)


The gene has 2 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Prrg4-201	ENSMUST00000028593.10	2872	226aa	Protein coding	CCDS16493	Q4KL73 Q8BGN6	TSL:1 GENCODE basic APPRIS is a system to annotate alternatively spliced transcripts based on a range of computational methods to identify the most functionally important transcript(s) of a gene. APPRIS P1
Prrg4-202	ENSMUST000000126824.1	751	145aa	Protein coding	-	A2A6Z6	CDS 3' incomplete TSL:2


The strategy is based on the design of *Prrg4-201* transcript,the transcription is shown below:



Genomic location distribution



Protein domain



If you have any questions, you are welcome to inquire.

Tel: 400-9660890

