

Mgat4d Cas9-KO Strategy

Designer: Yanhua Shen

Reviewer: Xueting Zhang

Design Date: 2020-2-14

Project Overview

Project Name

Mgat4d

Project type

Cas9-KO

Strain background

C57BL/6JGpt

Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Mgat4d* gene. The schematic diagram is as follows:



- The *Mgat4d* gene has 2 transcripts. According to the structure of *Mgat4d* gene, exon2-exon4 of *Mgat4d-201* (ENSMUST00000038692.5) transcript is recommended as the knockout region. The region contains 431bp coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Mgat4d* gene. The brief process is as follows: CRISPR/Cas9 system

- The *Mgat4d* gene is located on the Chr8. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This Strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of the gene knockout on gene transcription, RNA splicing and protein translation cannot be predicted at the existing technology level.

Gene information (NCBI)

Mgat4d MGAT4 family, member C [*Mus musculus* (house mouse)]

Gene ID: 67555, updated on 12-Oct-2019

Summary

Official Symbol	Mgat4d provided by MGI
Official Full Name	MGAT4 family, member C provided by MGI
Primary source	MGI:MGI:1914805
See related	Ensembl:ENSMUSG00000035057
Gene type	protein coding
RefSeq status	VALIDATED
Organism	Mus musculus
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	GnT1IP; GnT1IP-L; 4933434I20Rik
Expression	Restricted expression toward testis adult (RPKM 6.4) See more
Orthologs	human all

Genomic context

Location: 8; 8 C2

Exon count: 12


See Mgat4d in [Genome Data Viewer](#)

Transcript information (Ensembl)


The gene has 2 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Mgat4d-201	ENSMUST00000038692.5	1293	373aa	Protein coding	CCDS22451	Q9D4R2	TSL:1 GENCODE basic APPRIS P1
Mgat4d-202	ENSMUST00000139692.1	1690	No protein	lncRNA	-	-	TSL:1


The strategy is based on the design of *Mgat4d-201* transcript,The transcription is shown below




Genomic location distribution



Protein domain



Mouse phenotype description(MGI)



Phenotypes affected by the gene are marked in blue. Data quoted from MGI database(<http://www.informatics.jax.org/>).

If you have any questions, you are welcome to inquire.

Tel: 400-9660890

