

Haus7 Cas9-KO Strategy

Designer:

Yanhua Shen

Reviewer:

Xueting Zhang

Design Date:

2020-4-17

Project Overview

Project Name

Haus7

Project type

Cas9-KO

Strain background

C57BL/6JGpt

Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Haus7* gene. The schematic diagram is as follows:



- The *Haus7* gene has 3 transcripts. According to the structure of *Haus7* gene, exon3-exon4 of *Haus7-201* (ENSMUST00000033737.14) transcript is recommended as the knockout region. The region contains 130bp coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Haus7* gene. The brief process is as follows: CRISPR/Cas9 system

- Transcript 203 may not be affected.
- The *Haus7* gene is located on the ChrX. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This Strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of the gene knockout on gene transcription, RNA splicing and protein translation cannot be predicted at the existing technology level.

Gene information (NCBI)

Haus7 HAUS augmin-like complex, subunit 7 [Mus musculus (house mouse)]

Gene ID: 73738, updated on 13-Mar-2020

Summary




Official Symbol	Haus7 provided by MGI
Official Full Name	HAUS augmin-like complex, subunit 7 provided by MGI
Primary source	MGI:MGI:1920988
See related	Ensembl:ENSMUSG000000031371
Gene type	protein coding
RefSeq status	VALIDATED
Organism	Mus musculus
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	1110020L19Rik, Uchl5ip, Uip1
Expression	Ubiquitous expression in liver E14 (RPKM 18.5), CNS E11.5 (RPKM 17.8) and 28 other tissues See more
Orthologs	human all

Transcript information (Ensembl)


The gene has 3 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Haus7-201	ENSMUST00000033737.14	1380	364aa	Protein coding	CCDS30202	Q8BKT8	TSL:1 GENCODE basic APPRIS is a system to annotate alternatively spliced transcripts based on a range of computational methods to identify the most functionally important transcript(s) of a gene. APPRIS P4
Haus7-202	ENSMUST00000077243.4	1032	248aa	Protein coding	CCDS30203	Q8BKT8	TSL:1 GENCODE basic APPRIS is a system to annotate alternatively spliced transcripts based on a range of computational methods to identify the most functionally important transcript(s) of a gene. APPRIS ALT2
Haus7-203	ENSMUST00000131140.1	502	No protein	Processed transcript	-	-	TSL:3


The strategy is based on the design of *Haus7-201* transcript,The transcription is shown below



Genomic location distribution



Protein domain



If you have any questions, you are welcome to inquire.

Tel: 400-9660890

