

# Tc2n Cas9-KO Strategy

**Designer:** 

Daohua Xu

**Reviewer:** 

**Huimin Su** 

**Design Date:** 

2019-9-3

# **Project Overview**



**Project Name** 

Tc2n

**Project type** 

Cas9-KO


Strain background

C57BL/6JGpt

# **Knockout strategy**



This model will use CRISPR/Cas9 technology to edit the *Tc2n* gene. The schematic diagram is as follows:



### **Technical routes**



- ➤ The *Tc2n* gene has 4 transcripts. According to the structure of *Tc2n* gene, exon3-exon5 of *Tc2n-204*(ENSMUST00000162735.7) transcript is recommended as the knockout region. The region contains 491bp coding sequence.

  Knock out the region will result in disruption of protein function.
- > In this project we use CRISPR/Cas9 technology to modify Tc2n gene. The brief process is as follows: CRISPR/Cas9 system v

### **Notice**



- $\gt$  The Tc2n gene is located on the Chr12. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This Strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of the gene knockout on gene transcription, RNA splicing and protein translation cannot be predicted at the existing technology level.

## Gene information (NCBI)



#### Tc2n tandem C2 domains, nuclear [Mus musculus (house mouse)]

Gene ID: 74413, updated on 31-Jan-2019

#### Summary

☆ ?

Official Symbol Tc2n provided by MGI

Official Full Name tandem C2 domains, nuclear provided by MGI

Primary source MGI:MGI:1921663

See related Ensembl:ENSMUSG00000021187

Gene type protein coding
RefSeq status VALIDATED
Organism Mus musculus

Lineage Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha;

Muroidea; Muridae; Murinae; Mus; Mus

Also known as 4933406D09Rik, Mtac2d1, Tac2-N

Expression Biased expression in bladder adult (RPKM 2.4), testis adult (RPKM 0.5) and 6 other tissuesSee more

Orthologs <u>human</u> all


# Transcript information (Ensembl)



The gene has 4 transcripts, all transcripts are shown below:


Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Tc2n-204	ENSMUST00000162735.7	5290	489aa	Protein coding	CCDS36524	Q91XT6	TSL:1 GENCODE basic APPRIS P1
Tc2n-203	ENSMUST00000160830.1	2005	489aa	Protein coding	CCDS36524	Q91XT6	TSL:1 GENCODE basic APPRIS P1
Tc2n-202	ENSMUST00000160715.7	1386	<u>425aa</u>	Protein coding	-	E0CXQ6	TSL:5 GENCODE basic
Tc2n-201	ENSMUST00000110047.8	2587	489aa	Nonsense mediated decay	CCDS36524	Q91XT6	TSL:2

The strategy is based on the design of *Tc2n-204* transcript, The transcription is shown below




### Genomic location distribution





### Protein domain







If you have any questions, you are welcome to inquire. Tel: 400-9660890





