

Kcnj12 Cas9-KO Strategy

Designer: Huimin Su

Reviewer: Ruirui Zhang

Design Date: 2020/2/12

Project Overview



Project Name

Kcnj12

Project type

Cas9-KO


Strain background

C57BL/6JGpt

Knockout strategy



This model will use CRISPR/Cas9 technology to edit the Kcnj12 gene. The schematic diagram is as follows:



Technical routes



- ➤ The *Kcnj12* gene has 3 transcripts. According to the structure of *Kcnj12* gene, exon3 of *Kcnj12-201* (ENSMUST00000041944.8) transcript is recommended as the knockout region. The region contains most of the coding sequence. Knock out the region will result in disruption of protein function.
- ➤ In this project we use CRISPR/Cas9 technology to modify *Kcnj12* gene. The brief process is as follows: CRISPR/Cas9 system

Notice



- > According to the existing MGI data, Homozygotes for a targeted null mutation are viable and fertile with no detected abnormalities.
- The *Kcnj12* gene is located on the Chr11. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of the gene knockout on gene transcription, RNA splicing and protein translation cannot be predicted at the existing technology level.

Gene information (NCBI)



Kcnj12 potassium inwardly-rectifying channel, subfamily J, member 12 [Mus musculus (house mouse)]

Gene ID: 16515, updated on 12-Aug-2019

Summary

☆ ?

Official Symbol Kcnj12 provided by MGI

Official Full Name potassium inwardly-rectifying channel, subfamily J, member 12 provided by MGI

Primary source MGI:MGI:108495

See related Ensembl: ENSMUSG00000042529

Gene type protein coding
RefSeq status VALIDATED
Organism Mus musculus


Lineage Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia;

Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus

Also known as IRK2; IRK-2; Kir2.2; MB-IRK2

Expression Biased expression in cerebellum adult (RPKM 14.1), heart adult (RPKM 8.7) and 11 other tissues See more

Orthologs human all




Transcript information (Ensembl)



The gene has 3 transcripts, all transcripts are shown below:


Name 🍦	Transcript ID 🍦	bp 👙	Protein 4	Biotype	CCDS .	UniProt 4		Flags	4
Kcnj12-203	ENSMUST00000108717.2	2463	<u>427aa</u>	Protein coding	CCDS24805 €	<u>P52187</u> ₽	TSL:1	GENCODE basic	APPRIS P3
Kcnj12-201	ENSMUST00000041944.8	2459	<u>529aa</u>	Protein coding	CCDS59563 ₪	B1ATI1₽	TSL:1	GENCODE basic	APPRIS ALT2
Kcnj12-202	ENSMUST00000089184.10	2291	<u>427aa</u>	Protein coding	CCDS24805 €	<u>P52187</u> ₽	TSL:1	GENCODE basic	APPRIS P3

The strategy is based on the design of *Kcnj12-201* transcript, The transcription is shown below




Genomic location distribution





Protein domain





Mouse phenotype description(MGI)



Phenotypes affected by the gene are marked in blue. Data quoted from MGI database (http://www.informatics.jax.org/).

Homozygotes for a targeted null mutation are viable and fertile with no detected abnormalities.



If you have any questions, you are welcome to inquire. Tel: 400-9660890





