

***Sema4b* Cas9-KO Strategy**

Designer: Xueting Zhang

Reviewer: Yanhua Shen

Date: 2020-02-12

Project Overview

Project Name

Sema4b

Project type

Cas9-KO

Strain background

C57BL/6JGpt

Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Sema4b* gene. The schematic diagram is as follows:



- The *Sema4b* gene has 4 transcripts. According to the structure of *Sema4b* gene, exon3-exon13 of *Sema4b-201* (ENSMUST00000032754.8) transcript is recommended as the knockout region. The region contains 1531bp coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Sema4b* gene. The brief process is as follows: CRISPR/Cas9 system

- According to the existing MGI data, Mice homozygous for a gene trap allele exhibit normal cerebellar morphology. Mice homozygous for a knock-out allele exhibit enhanced memory response by way of increased IgE and IgG1 serum levels.
- The knockout region is near to the C-terminal of *Gm45206* gene, this strategy may influence the regulatory function of the C-terminal of *Gm45206* gene.
- The *Sema4b* gene is located on the Chr7. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This Strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of the gene knockout on gene transcription, RNA splicing and protein translation cannot be predicted at the existing technology level.

Gene information (NCBI)

Sema4b sema domain, immunoglobulin domain (Ig), transmembrane domain (TM) and short cytoplasmic domain, (semaphorin) 4B [*Mus musculus* (house mouse)]

Gene ID: 20352, updated on 6-Dec-2019

Summary

- Official Symbol** Sema4b provided by MGI
- Official Full Name** sema domain, immunoglobulin domain (Ig), transmembrane domain (TM) and short cytoplasmic domain, (semaphorin) 4B provided by MGI
- Primary source** MGI:MGI:107559
- See related** Ensembl:ENSMUSG00000030539
- Gene type** protein coding
- RefSeq status** VALIDATED
- Organism** [Mus musculus](#)
- Lineage** Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
- Also known as** SemC; Semac; sema C; mKIAA1745
- Expression** Broad expression in spleen adult (RPKM 31.2), duodenum adult (RPKM 28.9) and 27 other tissues [See more](#)
- Orthologs** [human](#) [all](#)

Genomic context

Location: 7; 7 D2 See Sema4b in [Genome Data Viewer](#)

Exon count: 16

Annotation release	Status	Assembly	Chr	Location
108	current	GRCm38.p6 (GCF_000001635.26)	7	NC_000073.6 (80186841..80226646)
Build 37.2	previous assembly	MGSCv37 (GCF_000001635.18)	7	NC_000073.5 (87331727..87371410)

Transcript information (Ensembl)


The gene has 4 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Sema4b-201	ENSMUST00000032754.8	3949	823aa	Protein coding	CCDS21391	Q62179	TSL:1 GENCODE basic APPRIS P1
Sema4b-204	ENSMUST00000205822.1	3760	823aa	Protein coding	CCDS21391	Q62179	TSL:1 GENCODE basic APPRIS P1
Sema4b-202	ENSMUST00000107383.7	2775	No protein	Retained intron	-	-	TSL:1
Sema4b-203	ENSMUST00000123023.2	1202	No protein	Retained intron	-	-	TSL:1

The strategy is based on the design of *Sema4b-201* transcript,The transcription is shown below




Genomic location distribution



Protein domain



Mouse phenotype description(MGI)



Phenotypes affected by the gene are marked in blue. Data quoted from MGI database(<http://www.informatics.jax.org/>).

According to the existing MGI data, Mice homozygous for a gene trap allele exhibit normal cerebellar morphology. Mice homozygous for a knock-out allele exhibit enhanced memory response by way of increased IgE and IgG1 serum levels.

If you have any questions, you are welcome to inquire.

Tel: 400-9660890

