

Gan Cas9-KO Strategy

Designer:

Daohua Xu

Reviewer:

Huimin Su

Design Date:

2020-2-25

Project Overview

Project Name

Gan

Project type


Cas9-KO

Strain background

C57BL/6JGpt

Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Gan* gene. The schematic diagram is as follows:



- The *Gan* gene has 2 transcripts. According to the structure of *Gan* gene, exon2-exon7 of *Gan-201* (ENSMUST00000064488.10) transcript is recommended as the knockout region. The region contains 1069bp coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Gan* gene. The brief process is as follows: CRISPR/Cas9 system w

- According to the existing MGI data, Null homozygotes display some muscular atrophy and motor neuron degeneration with the severity of these symptoms depending on genotype.
- The *Gan* gene is located on the Chr8. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This Strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of the gene knockout on gene transcription, RNA splicing and protein translation cannot be predicted at the existing technology level.

Gene information (NCBI)

Gan giant axonal neuropathy [Mus musculus (house mouse)]

Gene ID: 209239, updated on 31-Jan-2019

Summary




Official Symbol	Gan provided by MGI
Official Full Name	giant axonal neuropathy provided by MGI
Primary source	MGI:MGI:1890619
See related	Ensembl:ENSMUSG00000052557
Gene type	protein coding
RefSeq status	VALIDATED
Organism	Mus musculus
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	A330045G18, gigaxonin
Expression	Ubiquitous expression in lung adult (RPKM 4.0), adrenal adult (RPKM 3.8) and 28 other tissues See more
Orthologs	human all

Transcript information (Ensembl)


The gene has 2 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Gan-201	ENSMUST00000064488.10	2652	597aa	Protein coding	CCDS40490	Q8CA72	TSL:1 GENCODE basic APPRIS P1
Gan-202	ENSMUST00000162997.2	13461	597aa	Protein coding	-	F6TZU3	CDS 5' incomplete TSL:1


The strategy is based on the design of *Gan-201* transcript,The transcription is shown below




Genomic location distribution



Protein domain



Mouse phenotype description(MGI)



Phenotypes affected by the gene are marked in blue. Data quoted from MGI database(<http://www.informatics.jax.org/>).

According to the existing MGI data, Null homozygotes display some muscular atrophy and motor neuron degeneration with the severity of these symptoms depending on genotype.

If you have any questions, you are welcome to inquire.

Tel: 400-9660890

