

Asap2 Cas9-KO Strategy

Designer:Xueting Zhang

Reviwer: Yanhua Shen

Date:2020-02-20

Project Overview



Project Name

Asap2

Project type

Cas9-KO


Strain background

C57BL/6JGpt

Knockout strategy



This model will use CRISPR/Cas9 technology to edit the Asap2 gene. The schematic diagram is as follows:



Technical routes



- ➤ The Asap2 gene has 5 transcripts. According to the structure of Asap2 gene, exon2 of Asap2-202

 (ENSMUST00000064595.14) transcript is recommended as the knockout region. The region contains 73bp coding sequence.

 Knock out the region will result in disruption of protein function.
- ➤ In this project we use CRISPR/Cas9 technology to modify Asap2 gene. The brief process is as follows: CRISPR/Cas9 system

Notice



- ➤ The *Asap2* gene is located on the Chr12. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This Strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of the gene knockout on gene transcription, RNA splicing and protein translation cannot be predicted at the existing technology level.

Gene information (NCBI)



Asap2 ArfGAP with SH3 domain, ankyrin repeat and PH domain 2 [Mus musculus (house mouse)]

Gene ID: 211914, updated on 12-Aug-2019

Summary

☆ ?

Official Symbol Asap2 provided by MGI

Official Full Name ArfGAP with SH3 domain, ankyrin repeat and PH domain 2 provided by MGI

Primary source MGI:MGI:2685438

See related Ensembl: ENSMUSG00000052632

Gene type protein coding
RefSeq status VALIDATED
Organism Mus musculus

Lineage Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha;

Muroidea; Muridae; Murinae; Mus; Mus

Also known as PAP; PAG3; Ddef2; Gm592; Gm1523; 6530401G17Rik

Expression Ubiquitous expression in CNS E14 (RPKM 10.8), bladder adult (RPKM 10.1) and 27 other tissues See more

Genomic context



Location: 12; 12 A1.2-A1.3

See Asap2 in Genome Data Viewer

Exon count: 30

Annotation release	Status	Assembly	Chr	Location
108	current	GRCm38.p6 (GCF_000001635.26)	12	NC_000078.6 (2111136821270171)
Build 37.2	previous assembly	MGSCv37 (GCF_000001635.18)	12	NC_000078.5 (2111761721276032)


Transcript information (Ensembl)



The gene has 5 transcripts, all transcripts are shown below:


Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Asap2-202	ENSMUST00000064595.14	5678	1000aa	Protein coding	CCDS49037	E9PX52	TSL:5 GENCODE basic APPRIS P4
Asap2-204	ENSMUST00000101562.10	5552	<u>958aa</u>	Protein coding	CCDS56832	Q7SIG6	TSL:1 GENCODE basic APPRIS ALT2
Asap2-201	ENSMUST00000050990.9	3345	<u>955aa</u>	Protein coding	CCDS49036	D3YX85	TSL:5 GENCODE basic APPRIS ALT2
Asap2-203	ENSMUST00000090834.12	5097	809aa	Protein coding	-	Q7SIG6	TSL:5 GENCODE basic
Asap2-205	ENSMUST00000174530.1	1609	No protein	IncRNA		- 1	TSL:1

The strategy is based on the design of Asap2-202 transcript, The transcription is shown below




Genomic location distribution





Protein domain







If you have any questions, you are welcome to inquire. Tel: 400-9660890





