

Bag2 Cas9-KO Strategy

Designer: Xueting Zhang

Reviewer: Yanhua Shen

Date: 2020-02-24

Project Overview

Project Name

Bag2

Project type


Cas9-KO

Strain background

C57BL/6JGpt

Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Bag2* gene. The schematic diagram is as follows:



- The *Bag2* gene has 4 transcripts. According to the structure of *Bag2* gene, exon2-exon3 of *Bag2-201* (ENSMUST00000044691.8) transcript is recommended as the knockout region. The region contains most of the coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Bag2* gene. The brief process is as follows: CRISPR/Cas9 system v

Notice

- *Gm37905* gene will be deleted.
- Transcript *Bag2*-204 may not be affected.
- The *Bag2* gene is located on the Chr1. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This Strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of the gene knockout on gene transcription, RNA splicing and protein translation cannot be predicted at the existing technology level.

Gene information (NCBI)

Bag2 BCL2-associated athanogene 2 [*Mus musculus* (house mouse)]

Gene ID: 213539, updated on 12-Nov-2019

Summary

- Official Symbol** Bag2 provided by [MGI](#)
- Official Full Name** BCL2-associated athanogene 2 provided by [MGI](#)
- Primary source** [MGI:MGI:1891254](#)
- See related** [Ensembl:ENSMUSG00000042215](#)
- Gene type** protein coding
- RefSeq status** VALIDATED
- Organism** [Mus musculus](#)
- Lineage** Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
- Also known as** BC016230; 2610042A13Rik
- Expression** Ubiquitous expression in bladder adult (RPKM 26.9), placenta adult (RPKM 26.6) and 28 other tissues [See more](#)
- Orthologs** [human](#) [all](#)

Genomic context

Location: 1; 1 B

See Bag2 in [Genome Data Viewer](#)

Exon count: 4


Annotation release	Status	Assembly	Chr	Location
108	current	GRCm38.p6 (GCF_000001635.26)	1	NC_000067.6 (33745484..33757750, complement)
Build 37.2	previous assembly	MGSCv37 (GCF_000001635.18)	1	NC_000067.5 (33802329..33814595, complement)

Transcript information (Ensembl)


The gene has 4 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Bag2-201	ENSMUST00000044691.8	1860	210aa	Protein coding	CCDS14864	Q91YN9	TSL:1 GENCODE basic APPRIS P1
Bag2-203	ENSMUST00000187602.1	378	57aa	Protein coding	-	A0A087WNX9	TSL:2 GENCODE basic
Bag2-204	ENSMUST00000189741.1	955	No protein	Retained intron	-	-	TSL:NA
Bag2-202	ENSMUST00000155484.1	781	No protein	Retained intron	-	-	TSL:2


The strategy is based on the design of *Bag2-201* transcript,The transcription is shown below



Genomic location distribution



Protein domain



If you have any questions, you are welcome to inquire.

Tel: 400-9660890

